

Product data – L668

Product Name: LYPHAN boxes L668
pH 6,0 - 7,5

Description: LYPHAN boxes provide a quick measurement, when you require a precise measurement for liquids where you already know which general area of the pH scale they fall in. The test strips are impregnated with one indicator zone and comparison colors. These strips guarantee very high accuracy and precise readability.

Product description:	Product itself	Product with plastic box
Dimensions:	85 mm x 11mm	80 mm x 96 mm x 18 mm
Paper thickness:	0,27 - 0,30 mm	
Weight of the strip/box:	0,14 g	~ 65 g
Components:	filter paper; indicator; cotton dyes	
Indicator L669:	Bromoxylenol blue	

Indicator content per strip:
< 0,02 mg / strip (< 0,015 %)

Cotton dyes content per strip:
< 0,02 mg / strip (< 0,015 %)

Product picture:

