

Product data – LYPHAN boxes

Product Name: LYPHAN boxes L+ No. (e.g. L674; L666-8..)

Country of origin: Federal Republic of Germany

Customs tariff No.: 38220000

SKU information: 1 box = 200 strips

Product description: **Product itself** **Product with plastic box**

Dimensions: 85 mm x 11 mm 80 mm x 96 mm x 18 mm

Paper thickness: 0,27 - 0,30 mm

Weight of the strip/box: 0,14 g ~ 65 g with 200 strips

Components: filter paper; indicator; cotton dyes

No minimum order: 1 box with 200 strips

No UPC or EAN; no UN number; no UNSPSC; no barcode

no dangerous product; SDS attached

Storage temperature: room temperature

Shelf life: at least 12 month after receipt by consignee usable, if the strips stored under dry and dark conditions.

Product picture:

