

Product data – LYPHAN Rolls

Product Name: LYPHAN Roll R111/R3969/R4979/ R6999

Country of origin: Federal Republic of Germany

Customs tariff No.: 38220000

SKU information: 1 roll = 6 m

Product description: **Product itself** **Product with plastic box**

Length:	6 m	6 m
Diameter of the indicator role/box:	55 mm	max. 75 mm
Height:	10 mm	~ 18 mm
Paper thickness:	0,27 - 0,30 mm	
Weight of the roll/box:	~ 12 g	34,5 – 35 g

Components: filter paper; indicator
Indicator: indicator mixture; universal indicator
Indicator content per strip (~4 cm): < 0,05 mg (< 0,06 %)

No UPC or EAN; no UN number; no UNSPSC; no barcode

no dangerous product; SDS attached

Storage temperature: room temperature

Shelf life: at least 12 month after receipt by consignee usable, if the strips stored under dry and dark conditions.

Product picture:

