

TECHNICAL DATA SHEET

Part Number	049696
Description	Chest Freezer Racks



Benefits :

- ✓ Racks with drawers for all refrigerators and freezers.
- ✓ Stainless Steel
- ✓ 11 compartments for storage boxes.
- ✓ The freezer racks are equipped with colapsable handles on the sides.

Technical specifications :

- ✓ Weight : 1,54 kg
- ✓ Dimensions int. (Lxlxh) : 140 x 140 x 610 mm
- ✓ Rack Levels : 11 levels
- ✓ Box height up to: 54 mm
- ✓ Box dimensions up to : 136 x 136 mm
- ✓ Box capacity : 11

Sales unit :

1 Chest Freezer Racks
