

Artikel-Nr.: K25W
EAN: GLWAK25W
Manufacturer: GLW Storing Systems
GmbH

Color	white
Slides Capacity	25
Slide Width up to	1,3 mm
Autoclavable	no
Height	32 mm
Length	90 mm
Width	90 mm
Material	PS Polystyrene
Customs No.	3923 1090

GLW Slideboxes with 25 slots are made of polystyrene. The box accommodates temperature ranges from -80°C to +70°C. They are an optimal fit for standard microscope slides (76 mm x 26 mm x 1 mm). Non-standard microscope slides and round micro cover glasses of various sizes also fit perfectly in the box thanks to our accessory products (Item #: ZK20, ZK25 and ZK25W).

Each slidebox comes with slots numbered from 1 to 25 as well as a hinged cover lid. When closed, the boxes are stackable and can be secured in a stack by snapping them together. This enables the customer to save space when working or when storing the boxes. The hinged cover lid is detachable and can be removed and placed underneath the box. Due to this, the lid is not in the way and is not separated from the bottom half of the box.

They fit perfectly in the stainless steel working racks that we offer in our online shop.

The GLW Slideboxes are the ideal solution for storing, transporting and freezing slides.

