

Artikel-Nr.: K51

EAN: GLWAK51

Manufacturer: GLW Storing Systems GmbH

Color	transparent - autoclavable
Slides Capacity	50
Slide Width up to	1,3 mm
Autoclavable	yes
Height	31 mm
Length	172 mm
Width	83 mm
Material	PC Polycarbonate
Customs No.	3923 1090

GLW Slideboxes with 50 slots are available in two different materials, polystyrene and polycarbonate. The polystyrene boxes accommodate temperature ranges from -80°C to +70°C and the polycarbonate boxes accommodate temperature ranges from -150°C to +130°C and are **autoclavable**. They are an optimal fit for standard microscope slides (76 mm x 26 mm x 1 mm). Non-standard slides and cover glasses of various sizes also fit perfectly in the box thanks to our accessory products (Item #: ZK20, ZK25, ZK25W and ZKH301). For example, we offer Slidebox Height Extenders (Item #: ZKH301) in our online shop, which enable you to store larger microscope slides (up to 30 mm larger) in the box. Each slidebox comes with slots numbered from 1 to 50 as well as a snap-on lid. The box comes with a numbered index card.

When closed, the boxes are stackable and can be secured in a stack by snapping them together. This enables the customer to save space when working or when storing the boxes. The hinged cover lid is detachable and can be removed and placed underneath the box. Due to this, the lid is not in the way and is not separated from the bottom half of the box.

They fit perfectly in the stainless steel working racks that we offer in our online shop.

The GLW Slideboxes are the ideal solution for storing, transporting and freezing slides.

