


DACOS

Cryo Vials • Pipette Tips • Blood-Sep Filters

MAXXLINE 9X9 CRYO BOX DATASHEET



MAXXLINE 9x9 Cryo Boxes

The MAXXLINE 9x9 Cryogenic Boxes are a practical designed storage solution for cryogenic vials. With features like graduated alpha-numeric coding, smart closing identification and high-quality durable resin, these boxes are made with the purpose of fast and easy handling.

INTENDED USE:

General laboratory product, solely intended for general laboratory use by competent staff.

Item no.	CM81-[x]	CM81-[x]-4	CM81-[x]-5
FULL ITEM NO.			
N - Nature	CM81-N	CM81-N-4	CM81-N-5
B - Blue	CM81-B	CM81-B-4	CM81-B-5
G - Green	CM81-G	CM81-G-4	CM81-G-5
R - Red	CM81-R	CM81-R-4	CM81-R-5
Y - Yellow	CM81-Y	CM81-Y-4	CM81-Y-5
SPECIFICATIONS			
Vial size	1.0 + 2.0 ml	4.0 ml	5.0 ml
Height (mm)	50	80	90
Length x width (mm)	130x130		
Material	PP (Polypropylen)		
Temperature range	-90° to 121°C <i>Note: Internally tested to withstand temperatures down to -150°C.</i>		
Autoclavable	Yes		
FEATURES			
Lid	Clear lid - Enables perfect, undisturbed view of samples. Alpha-numerical system - Molded in for easy sample identification.		
Base	Colors - The base is available in 5 different colors. Ventilation - 4 openings for ventilation.		
PACKAGING & STORAGE			
Pcs/bag	105	90	60
LOT no.	LOT number printed on each product pouch and box label.		
Storage	Store in closed, original packaging in a dry place at room temperature (10-30°C). Do not expose to direct sunlight.		

+45 76307600 • info@dacos.dk • www.dacos.dk

DOMINIQUE DUTSCHER SAS

