

Appendix to EC DECLARATION OF CONFORMITY for Finnpiquette™ Digital, F1, F2 and F3

We hereby certify that Finnpiquette™ Pipettes (Finnpiquette Digital, F1, F2 and F3) and all Finntip™ Pipette Tips listed in Appendix 1, are marked with CE symbol and they comply with the European Union Directive 98/79/EC.

EC Declaration of Conformity for Finnpiquette Digital, F1, F2 and F3, with ID numbers listed below are valid for the products listed in this Appendix 1.

Finnpiquette	EC DoC
Finnpiquette Digital	15163820-07
Finnpiquette F1	15283850-05
Finnpiquette F2	15283840-05
Finnpiquette F3	15283820-05


The CE mark covers the combination of pipettes and tips and accessories defined in Appendix 1. Any other pipette – tip combination has to be validated by the user.

Our Quality System has been found to conform to Management System Standards ISO 9001:2008 and ISO 13485:2003.

Vantaa, June 29, 2016



 Jari Mäkilaakso
 Regulative Manager
 Liquid Handling Consumables
 Thermo Fisher Scientific Oy





 Marina Elf
 Quality and Compliance Manager
 Liquid Handling Consumables
 Thermo Fisher Scientific Oy

Finnpipette™ Pipettes

Catalogue number	Description
4500000	FP-Digital 0.2-2µl
4500010	FP-Digital 0.5-10µl µTip
4500020	FP-Digital 0.5 - 10µl
4500030	FP-Digital 5-40µl
4500050	FP-Digital 200-1000µl
4500060	FP-Digital 1-5ml
4500070	FP-Digital 2-10ml
4500080	FP-Digital 2-20µl
4500090	FP-Digital 20-200µl
4500100	FP-Digital 5-50µl
4500110	FP-Digital 10-100µl
4500120	FP-Digital 100-1000µl
4510000	FP-Digital MCP8 0.5-10µl
4510010	FP-Digital MCP12 0.5-10µl
4510020	FP-Digital MCP8 5-50µl
4510030	FP-Digital MCP8 50-300µl
4510040	FP-Digital MCP12 5-50µl
4510050	FP-Digital MCP12 50-300µl
4510070	FP-Digital MCP16 5-50µl
4640000	Finnpipette F3 1-10µl µTip
4640010	Finnpipette F3 1-10µl
4640020	Finnpipette F3 2-20µl µTip
4640030	Finnpipette F3 2-20µl
4640040	Finnpipette F3 10-100µl
4640050	Finnpipette F3 20-200µl
4640060	Finnpipette F3 100-1000µl
4640070	Finnpipette F3 1-10ml
4640080	Finnpipette F3 5-50µl µTip
4640090	Finnpipette F3 5-50µl
4640100	Finnpipette F3 0.5-5ml
4640110	Finnpipette F3 0.2-2µl µTip
4641010	Finnpipette F1 0.2-2µl
4641010N	Finnpipette F1 0.2-2µl
4641020	Finnpipette F1 0.5-5µl
4641020N	Finnpipette F1 0.5-5µl
4641030	Finnpipette F1 1-10µl µTip
4641030N	Finnpipette F1 1-10µl µTip
4641040	Finnpipette F1 1-10µl
4641040N	Finnpipette F1 1-10µl
4641050	Finnpipette F1 2-20µl µTip
4641050N	Finnpipette F1 2-20µl µTip
4641060	Finnpipette F1 2-20µl
4641060N	Finnpipette F1 2-20µl
4641070	Finnpipette F1 10-100µl
4641070N	Finnpipette F1 10-100µl
4641080	Finnpipette F1 20-200µl
4641080N	Finnpipette F1 20-200µl
4641090	Finnpipette F1 30-300µl
4641090N	Finnpipette F1 30-300µl

4641100	Finnpipette F1 100-1000µl
4641100N	Finnpipette F1 100-1000µl
4641110	Finnpipette F1 0.5-5ml
4641110N	Finnpipette F1 0.5-5ml
4641120	Finnpipette F1 1-10ml
4641120N	Finnpipette F1 1-10ml
4641130	Finnpipette F1 5-50µl µTip
4641130N	Finnpipette F1 5-50µl µTip
4641140	Finnpipette F1 5-50µl
4641140N	Finnpipette F1 5-50µl
4642010	Finnpipette F2 0.2-2µl
4642020	Finnpipette F2 0.5-5µl
4642030	Finnpipette F2 1-10µl µTip
4642040	Finnpipette F2 1-10µl
4642050	Finnpipette F2 2-20µl µTip
4642060	Finnpipette F2 2-20µl
4642070	Finnpipette F2 10-100µl
4642080	Finnpipette F2 20-200µl
4642090	Finnpipette F2 100-1000µl
4642100	Finnpipette F2 0.5-5ml
4642110	Finnpipette F2 1-10ml
4642120	Finnpipette F2 5-50µl µTip
4642130	Finnpipette F2 5-50µl
4650000	Finnpipette F3 Fixed 1µl µTip
4650010	Finnpipette F3 Fixed 5µl µTip
4650020	Finnpipette F3 Fixed 10µl
4650030	Finnpipette F3 Fixed 25µl
4650040	Finnpipette F3 Fixed 50µl
4650050	Finnpipette F3 Fixed 100µl
4650060	Finnpipette F3 Fixed 250µl
4650070	Finnpipette F3 Fixed 500µl
4650080	Finnpipette F3 Fixed 1000µl
4650090	Finnpipette F3 Fixed 20µl
4650100	Finnpipette F3 Fixed 200µl
4651000	Finnpipette F1 Fixed 1µl µTip
4651010	Finnpipette F1 Fixed 5µl µTip
4651020	Finnpipette F1 Fixed 10µl
4651020N	Finnpipette F1 Fixed 10µl
4651030	Finnpipette F1 Fixed 25µl
4651030N	Finnpipette F1 Fixed 25µl
4651040	Finnpipette F1 Fixed 50µl
4651040N	Finnpipette F1 Fixed 50µl
4651050	Finnpipette F1 Fixed 100µl
4651050N	Finnpipette F1 Fixed 100µl
4651060	Finnpipette F1 Fixed 250µl
4651060N	Finnpipette F1 Fixed 250µl
4651070	Finnpipette F1 Fixed 500µl
4651070N	Finnpipette F1 Fixed 500µl
4651080	Finnpipette F1 Fixed 1000µl
4651080N	Finnpipette F1 Fixed 1000µl
4651090	Finnpipette F1 Fixed 2ml
4651090N	Finnpipette F1 Fixed 2ml

4651100	Finnpipette F1 Fixed 3ml
4651100N	Finnpipette F1 Fixed 3ml
4651110	Finnpipette F1 Fixed 5ml
4651110N	Finnpipette F1 Fixed 5ml
4651120	Finnpipette F1 Fixed 10ml
4651120N	Finnpipette F1 Fixed 10ml
4651130	Finnpipette F1 Fixed 20µl
4651130N	Finnpipette F1 Fixed 20µl
4651140	Finnpipette F1 Fixed 200µl
4651140N	Finnpipette F1 Fixed 200µl
4652000	Finnpipette F2 Fixed 1µl µTip
4652010	Finnpipette F2 Fixed 5µl µTip
4652020	Finnpipette F2 Fixed 10µl
4652030	Finnpipette F2 Fixed 25µl
4652040	Finnpipette F2 Fixed 50µl
4652050	Finnpipette F2 Fixed 100µl
4652060	Finnpipette F2 Fixed 250µl
4652070	Finnpipette F2 Fixed 500µl
4652080	Finnpipette F2 Fixed 1000µl
4652090	Finnpipette F2 Fixed 2ml
4652100	Finnpipette F2 Fixed 3ml
4652110	Finnpipette F2 Fixed 5ml
4652120	Finnpipette F2 Fixed 10ml
4652130	Finnpipette F2 Fixed 20µl
4652140	Finnpipette F2 Fixed 200µl
4652150	Finnpipette F2 Fixed 4ml
4660010	Finnpipette F3 MCP8 5-50µl
4660020	Finnpipette F3 MCP8 30-300µl
4660040	Finnpipette F3 MCP12 5-50µl
4660050	Finnpipette F3 MCP12 30-300µl
4661000	Finnpipette F1 MCP8 1-10µl
4661000N	Finnpipette F1 MCP8 1-10µl
4661010	Finnpipette F1 MCP8 5-50µl
4661010N	Finnpipette F1 MCP8 5-50µl
4661020	Finnpipette F1 MCP8 10-100µl
4661020N	Finnpipette F1 MCP8 10-100µl
4661030	Finnpipette F1 MCP8 30-300µl
4661030N	Finnpipette F1 MCP8 30-300µl
4661040	Finnpipette F1 MCP12 1-10µl
4661040N	Finnpipette F1 MCP12 1-10µl
4661050	Finnpipette F1 MCP12 5-50µl
4661050N	Finnpipette F1 MCP12 5-50µl
4661060	Finnpipette F1 MCP12 10-100µl
4661060N	Finnpipette F1 MCP12 10-100µl
4661070	Finnpipette F1 MCP12 30-300µl
4661070N	Finnpipette F1 MCP12 30-300µl
4661080	Finnpipette F1 MCP16 1-10µl
4661080N	Finnpipette F1 MCP16 1-10µl
4661090	Finnpipette F1 MCP16 5-50µl
4661090N	Finnpipette F1 MCP16 5-50µl
4662000	Finnpipette F2 MCP8 1-10µl
4662010	Finnpipette F2 MCP8 5-50µl

4662020	Finnpipette F2 MCP8 10-100µl
4662030	Finnpipette F2 MCP8 30-300µl
4662040	Finnpipette F2 MCP12 1-10µl
4662050	Finnpipette F2 MCP12 5-50µl
4662060	Finnpipette F2 MCP12 10-100µl
4662070	Finnpipette F2 MCP12 30-300µl
4662080	Finnpipette F2 MCP16 1-10µl
4662090	Finnpipette F2 MCP16 5-50µl

Finntip™ Pipette tips

Catalogue number	Description
9400100	Finntip 200 Ext, 400/box
9400130	Finntip 200 Ext, 10x96/rack
9400133	Finntip 200 Ext, 10x96/rack, sterile
9400220	Finntip 250 Universal, 20000/bulk
9400230	Finntip 250 Universal, 1000/bag
9400250	Finntip 250 Universal, 500/box
9400260	Finntip 250 Universal, 10x96/rack
9400263	Finntip 250 Universal, 10x96/rack, sterile
9400266	Finntip 250 Universal, refill starter kit, 10x96/refill+4x96/rack
9400267	Finntip 250 Universal, 20x96/refill kit
9400300	Finntip 10, 10x96/rack
9400303	Finntip 10, 10x96/rack, sterile
9400310	Finntip 10, 1000/bag
9400326	Finntip 10, refill starter kit, 10x192/refill+4x192/rack
9400327	Finntip 10, 20x192/refill kit
9400360	Finntip 50, 1000/bag
9400370	Finntip 50, 10x384/rack
9400373	Finntip 50, 10x384/rack, sterile
9400610	Finntip 20, 10x384/rack
9400613	Finntip 20, 10x384/rack, sterile
9400620	Finntip 20, 1000/bag
9401010	Finntip 1000, 10000/bulk
9401030	Finntip 1000, 1000/bag
9401070	Finntip 1000, 200/box
9401110	Finntip 1000, 10x96/rack
9401113	Finntip 1000, 10x96/rack, sterile
9401115	Finntip 1000, 5x192/refill pack
9401240	Finntip 300, 20000/bulk
9401250	Finntip 300, 10x96/rack
9401253	Finntip 300, 10x96/rack, sterile
9401255	Finntip 300, 10x192/refill pack
9401260	Finntip 300, 1000/bag
9401410	Finntip 1000 Ext, 1000/bag
9401420	Finntip 1000 Ext, 5x96/rack
9401423	Finntip 1000 Ext, 5x96/rack, sterile
9402010	Finntip 5 ml, 3000/bulk
9402030	Finntip 5 ml, 500/bag
9402050	Finntip 5 ml, 75/box
9402070	Finntip 5 ml, 5x54/rack

9402073	Finntip 5 ml, 5x54/rack, sterile
9402150	Finntip 10 ml, 40/box
9402151	Finntip 10 ml, 100/bag
9402160	Finntip 10 ml, 5x24/rack
9402163	Finntip 10 ml, 5x24/rack, sterile
9405020	Finntip 250 Wide, 1000/bag
9405030	Finntip 250 Wide, 22000/bulk
9405050	Finntip 1000 Wide, 400/bag
9405060	Finntip 1000 Wide, 6500/bulk
9405120	Finntip 250 Wide, 10x96/rack
9405123	Finntip 250 Wide, 10x96/rack, sterile
9405160	Finntip 1000 Wide, 10x96/rack
9405163	Finntip 1000 Wide, 10x96/rack, sterile
94052000	Finntip Filter 10µl Micro, 10x96/rack, sterile
94052020	Finntip Filter 20µl, 10x384/rack, sterile
94052060	Finntip Filter 50µl, 10x384/rack, sterile
94052100	Finntip Filter 10µl, 10x96/rack, sterile
94052150	Finntip Filter 20µl, 10x96/rack, sterile
94052160	Finntip Filter 30µl, 10x96/rack, sterile
94052200	Finntip Filter 100µl, 10x96/rack, sterile
94052300	Finntip Filter 200µl, 10x96/rack, sterile
94052310	Finntip Filter 100µl Ext, 10x96/rack, sterile
94052320	Finntip Filter 200µl Ext, 10x96/rack, sterile
94052350	Finntip Filter 300µl, 10x96/rack, sterile
94052410	Finntip Filter 1000µl, 10x96/rack, sterile
94052430	Finntip Filter 1000µl Ext, 5x96/rack, sterile
94052450	Finntip Filter 1500µl 10x96/rack, sterile
94052550	Finntip Filter 5ml, 5x54/rack, sterile
94052600	Finntip Filter 10ml, 5x24/rack, sterile
94053100	Finntip BioCon 10, 100/box, sterile
94053200	Finntip BioCon 250, 100/box, sterile
94053300	Finntip BioCon 1000, 50/box, sterile
94056380	Finntip Flex Filter 200µl, 10x96/rack, sterile
94056510	Finntip Flex Filter 30µl, 10x96/rack, sterile
94056520	Finntip Flex Filter 100µl, 10x96/rack, sterile
94056580	Finntip Flex Filter 300µl, 10x96/rack, sterile
94056710	Finntip Flex Filter 1000µl, 10x96/rack, sterile
94056810	Finntip Flex Filter 1200µl, 10x96/rack, sterile
94056970	Finntip Flex Filter 10ml Ext, 50/bag, sterile
94056980	Finntip Flex Filter 10µl, 10x96/rack, sterile
94060100	Finntip Flex 10, 10x96/rack
94060103	Finntip Flex 10, 10x96/rack, sterile
94060116	Finntip Flex 10, refill starter kit, 20x96/refill+1x96/rack
94060117	Finntip Flex 10, 20x96/refill kit
94060120	Finntip Flex 10, 1000/bag
94060310	Finntip Flex 200, 10x96/rack
94060313	Finntip Flex 200, 10x96/rack, sterile
94060316	Finntip Flex 200, refill starter kit, 20x96/refill+1x96/rack
94060317	Finntip Flex 200, 20x96/refill kit
94060320	Finntip Flex 200, 1000/bag
94060510	Finntip Flex 300, 10x96/rack
94060513	Finntip Flex 300, 10x96/rack, sterile



94060516	Finntip Flex 300, refill starter kit, 20x96/refill+1x96/rack
94060517	Finntip Flex 300, 20x96/refill kit
94060520	Finntip Flex 300, 1000/bag
94060710	Finntip Flex 1000, 10x96/rack
94060713	Finntip Flex 1000, 10x96/rack, sterile
94060716	Finntip Flex 1000, refill starter kit, 16x96/refill+1x96/rack
94060717	Finntip Flex 1000, 16x96/refill kit
94060720	Finntip Flex 1000, 1000/bag
94060810	Finntip Flex 1200, 10x96/rack
94060813	Finntip Flex 1200, 10x96/rack, sterile
94060816	Finntip Flex 1200, refill starter kit, 16x96/refill+1x96/rack
94060817	Finntip Flex 1200, 16x96/refill kit
94060820	Finntip Flex 1200, 1000/bag
94060970	Finntip Flex 10 ml Ext, 100/bag
94060973	Finntip Flex 10 ml Ext, 50/bag, sterile

Finnpipette™ Accessories

9420290	Pipette stand, gray, for 6 pipettes
9420400	F-stand, Linear Stand, white, 6 positions
9420320	Pipette Mini-stand, gray, for 3 pipettes
9420390	Multichannel Pipette Stand, 1 position
9420030	Multichannel Pipette Stand, white

