

Appendix to Declaration of Conformity: Finnpipette Digital, F1, F2 and F3

We hereby certify that Finnpipettes (Finnpipette Digital, F1, F2 and F3) and all Finntips listed in Appendix 1, are marked with CE symbol and they comply with the European Union Directive 98/79/EC.

EC Declarations of Conformity for Finnpipette Digital, F1, F2 and F3, with ID Numbers listed below and signed Oct 18 2012, are valid for the products listed in Appendix 1 and Appendix 2.

<u>Finnpipette</u>	<u>EC DoC</u>
Finnpipette Digital	15163820-04
Finnpipette F1	15283850-02
Finnpipette F2	15283840-02
Finnpipette F3	15283820-04

The CE mark covers the combination of above mentioned pipettes, tips defined in Appendix 1 and accessories defined in Appendix 2. Any other pipette – tip combination has to be validated by the user.

Our Quality System has been found to conform to the Management System Standards ISO 9001:2008 and ISO 13485:2003.

In Vantaa, on 23rd of October 2014



Marina Elf
Quality Manager
Liquid Handling Consumables
Thermo Fisher Scientific Oy

APPENDIX 1

Finntips

Catalogue number	Description
9400100	Finntip 200 Ext, 400/box
9400130	Finntip 200 Ext, 10x96/rack
9400133	Finntip 200 Ext, 10x96/rack, sterile
9400220	Finntip 250 Universal, 20000/bulk
9400230	Finntip 250 Universal, 1000/bag
9400250	Finntip 250 Universal, 500/box
9400260	Finntip 250 Universal, 10x96/rack
9400263	Finntip 250 Universal, 10x96/rack, sterile
9400266	Finntip 250 Universal, refill starter kit, 10x96/refill+4x96/rack
9400267	Finntip 250 Universal, 20x96/refill kit
9400300	Finntip 10, 10x96/rack
9400303	Finntip 10, 10x96/rack, sterile
9400310	Finntip 10, 1000/bag
9400326	Finntip 10, refill starter kit, 10x192/refill+4x192/rack
9400327	Finntip 10, 20x192/refill kit
9400360	Finntip 50, 1000/bag
9400370	Finntip 50, 10x384/rack
9400373	Finntip 50, 10x384/rack, sterile
9400610	Finntip 20, 10x384/rack
9400613	Finntip 20, 10x384/rack, sterile
9400620	Finntip 20, 1000/bag
9401010	Finntip 1000, 10000/bulk
9401030	Finntip 1000, 1000/bag
9401070	Finntip 1000, 200/box
9401110	Finntip 1000, 10x96/rack
9401113	Finntip 1000, 10x96/rack, sterile
9401115	Finntip 1000, 5x192/refill pack
9401240	Finntip 300, 20000/bulk
9401250	Finntip 300, 10x96/rack
9401253	Finntip 300, 10x96/rack, sterile
9401255	Finntip 300, 10x192/refill pack
9401260	Finntip 300, 1000/bag
9401410	Finntip 1000 Ext, 1000/bag
9401420	Finntip 1000 Ext, 5x96/rack
9401423	Finntip 1000 Ext, 5x96/rack, sterile
9402010	Finntip 5 ml, 3000/bulk
9402030	Finntip 5 ml, 500/bag
9402050	Finntip 5 ml, 75/box
9402070	Finntip 5 ml, 5x54/rack
9402073	Finntip 5 ml, 5x54/rack, sterile
9402150	Finntip 10 ml, 40/box
9402151	Finntip 10 ml, 100/bag
9402160	Finntip 10 ml, 5x24/rack
9402163	Finntip 10 ml, 5x24/rack, sterile
9405020	Finntip 250 Wide, 1000/bag

9405030	Finntip 250 Wide, 22000/bulk
9405050	Finntip 1000 Wide, 400/bag
9405060	Finntip 1000 Wide, 6500/bulk
9405120	Finntip 250 Wide, 10x96/rack
9405123	Finntip 250 Wide, 10x96/rack, sterile
9405160	Finntip 1000 Wide, 10x96/rack
9405163	Finntip 1000 Wide, 10x96/rack, sterile
94052000	Finntip Filter 10µl Micro, 10x96/rack, sterile
94052020	Finntip Filter 20µl, 10x384/rack, sterile
94052060	Finntip Filter 50µl, 10x384/rack, sterile
94052100	Finntip Filter 10µl, 10x96/rack, sterile
94052150	Finntip Filter 20µl, 10x96/rack, sterile
94052160	Finntip Filter 30µl, 10x96/rack, sterile
94052200	Finntip Filter 100µl, 10x96/rack, sterile
94052300	Finntip Filter 200µl, 10x96/rack, sterile
94052310	Finntip Filter 100µl Ext, 10x96/rack, sterile
94052320	Finntip Filter 200µl Ext, 10x96/rack, sterile
94052350	Finntip Filter 300µl, 10x96/rack, sterile
94052410	Finntip Filter 1000µl, 10x96/rack, sterile
94052430	Finntip Filter 1000µl Ext, 5x96/rack, sterile
94052450	Finntip Filter 1500µl 10x96/rack, sterile
94052550	Finntip Filter 5ml, 5x54/rack, sterile
94052600	Finntip Filter 10ml, 5x24/rack, sterile
94053100	Finntip BioCon 10, 100/box, sterile
94053200	Finntip BioCon 250, 100/box, sterile
94053300	Finntip BioCon 1000, 50/box, sterile
94056380	Finntip Flex Filter 200µl, 10x96/rack, sterile
94056510	Finntip Flex Filter 30µl, 10x96/rack, sterile
94056520	Finntip Flex Filter 100µl, 10x96/rack, sterile
94056580	Finntip Flex Filter 300µl, 10x96/rack, sterile
94056710	Finntip Flex Filter 1000µl, 10x96/rack, sterile
94056810	Finntip Flex Filter 1200µl, 10x96/rack, sterile
94056970	Finntip Flex Filter 10ml Ext, 50/bag, sterile
94056980	Finntip Flex Filter 10µl, 10x96/rack, sterile
94060100	Finntip Flex 10, 10x96/rack
94060103	Finntip Flex 10, 10x96/rack, sterile
94060116	Finntip Flex 10, refill starter kit, 20x96/refill+1x96/rack
94060117	Finntip Flex 10, 20x96/refill kit
94060120	Finntip Flex 10, 1000/bag
94060310	Finntip Flex 200, 10x96/rack
94060313	Finntip Flex 200, 10x96/rack, sterile
94060316	Finntip Flex 200, refill starter kit, 20x96/refill+1x96/rack
94060317	Finntip Flex 200, 20x96/refill kit
94060320	Finntip Flex 200, 1000/bag
94060510	Finntip Flex 300, 10x96/rack
94060513	Finntip Flex 300, 10x96/rack, sterile
94060516	Finntip Flex 300, refill starter kit, 20x96/refill+1x96/rack
94060517	Finntip Flex 300, 20x96/refill kit
94060520	Finntip Flex 300, 1000/bag
94060710	Finntip Flex 1000, 10x96/rack

94060713	Finntip Flex 1000, 10x96/rack, sterile
94060716	Finntip Flex 1000, refill starter kit, 16x96/refill+1x96/rack
94060717	Finntip Flex 1000, 16x96/refill kit
94060720	Finntip Flex 1000, 1000/bag
94060810	Finntip Flex 1200, 10x96/rack
94060813	Finntip Flex 1200, 10x96/rack, sterile
94060816	Finntip Flex 1200, refill starter kit, 16x96/refill+1x96/rack
94060817	Finntip Flex 1200, 16x96/refill kit
94060820	Finntip Flex 1200, 1000/bag
94060970	Finntip Flex 10 ml Ext, 100/bag
94060973	Finntip Flex 10 ml Ext, 50/bag, sterile

APPENDIX 2

Accessories

9420290	Pipette stand, gray, for 6 pipettes
9420400	F-stand, Linear Stand, white, 6 positions
9420320	Pipette Mini-stand, gray, for 3 pipettes
9420390	Multichannel Pipette Stand, 1 position
9420030	Multichannel Pipette Stand, white