

# Terumo Syringe

3-part Syringe



## Product Information

NEW PRODUCT CODES!



### Disposable syringes for single use

Every single detail is the result of Terumo's long time experience in this field which is appreciated by users throughout the world.

- Siliconisation for optimal plunger movement
- Firm and definite end-stop
- Clear scale markings to enable precise dosage
- Transparent barrel for visualisation
- Finger grips support easy handling

#### Sterilization

With e-beam (for ,SS'-codes) or ethylene oxide (for ,MDSS'-codes) to a Sterility Assurance Level of at least  $10^{-6}$

#### Storage

- Shelf life of 5 years
- The batch number and expiry date are embossed in the film of the blister pack
- Do not store at extreme temperature and humidity

Description	Raw Materials
1 Barrel	Polypropylene
2 Plunger	Polypropylene
3 Gasket	Thermoplastic elastomer (for ,SS'-codes) Polyisoprene rubber (for ,MDSS'-codes)
Lubricant	Silicone oil

Does not contain components made of PVC or natural rubber latex.

#### Precautions

- Do not use if package is damaged
- Do not reuse
- Dispose of sharps after single use according to facility protocol to avoid risk of infection

# Terumo Syringe

3-part Syringe



## Your syringes of choice in different sizes and tip types

**NEW PRODUCT CODES!**


### Product Range

Product Code	Nominal Capacity	Tip	Replacement for
MDSS01SE	1 ml	Concentric Luer Slip	SS+01T1
SS*02LE1	2 ml	Concentric Luer Lock	-
MDSS03SE	3 ml	Concentric Luer Slip	SS*02SE1, SS+03S
MDSS03LE	3 ml	Concentric Luer Lock	SS*03L, SS+03L1
MDSS05SE	5 ml	Concentric Luer Slip	SS*05SE1, SS+05S
SS*05LE1	5 ml	Concentric Luer Lock	-
MDSS10SE	10 ml	Eccentric Luer Slip	SS+10ES1
SS*10LE1	10 ml	Concentric Luer Lock	-
MDSS20ESE	20 ml	Eccentric Luer Slip	SS+20ES1
SS+20L1	20 ml	Concentric Luer Lock	-
SS+30L1MP	30 ml	Concentric Luer Lock	-
SS+30ES1	30 ml	Eccentric Luer Slip	-
SS+50L1	50 ml	Concentric Luer Lock	-
SS+50ES1	50 ml	Eccentric Luer Slip	-
SS+50C1	50 ml	Catheter Tip	-

### Packaging

Unit Box	Shipping Carton		
	Quantity	Gross Weight	Dimensions (L x W x H)
100	3000	16.60 kg	445 x 445 x 445 mm
100	2400	6.00 kg	578 x 386 x 436 mm
100	1800	10.90 kg	445 x 445 x 445 mm
100	1800	10.90 kg	445 x 445 x 445 mm
100	1800	14.90 kg	655 x 445 x 425 mm
100	1600	6.25 kg	544 x 404 x 444 mm
100	1200	14.60 kg	565 x 425 x 485 mm
100	800	15.00 kg	580 x 350 x 377 mm
50	400	14.60 kg	465 x 385 x 455 mm
50	600	10.15 kg	570 x 400 x 366 mm
25	200	6.32 kg	552 x 386 x 238 mm
25	200	6.20 kg	552 x 386 x 238 mm
25	100	5.10 kg	466 x 312 x 268 mm
25	100	5.10 kg	466 x 312 x 268 mm
25	100	5.40 kg	466 x 312 x 268 mm

Single unit: Blister pack (film/film)



### Regulatory Status

CE Certification – Class IIa medical device

Notified body TÜV Rheinland 0197

Manufacturer for .SS+/-codes TERUMO (Philippines) CORPORATION, 124 East Main Avenue, Laguna Technopark, Binan, Philippines  
for .SS\*-codes TERUMO CORPORATION, 2-44-1 Hatagaya, Shibuya-Ku, Toyko 151-0072, Japan  
for .MDSS'-codes Yangzhou Medline Industry Co., Ltd. (China)

EC representative for .SS'-codes TERUMO EUROPE N.V. – Leuven, Belgium  
for .MDSS'-codes Shanghai International Corp. GmbH (Europe) – Hamburg, Germany

Distributor for .MDSS'-codes TERUMO EUROPE N.V., Interleuvenlaan 40, 3001 Leuven, Belgium

Choose from a  
range of volumes  
from 1 ml to 50 ml

Syringes are available in all  
common tip types.



Terumo Corporation  
+81 3 6742 8500

Terumo Europe N.V.  
+32 16 38 12 11

Your Terumo partner

All brand names are trademarks or registered trademarks of TERUMO CORPORATION and their respective owners.  
Published by Terumo Europe N.V.

Version A  
Please note that some products may not be available in all regions, check with your local Terumo contact.