



Celebrating over
60 years of service

WORLDWIDE SERVICE TO EDUCATION • HEALTH • AGRICULTURE • INDUSTRY

901 JANESVILLE AVENUE FORT ATKINSON, WISCONSIN 53538-0901
PHONE 920-563-2446 FAX 920-563-8296
www.eNASCO.com E-MAIL info@eNASCO.com

**Whirl-Pak® Sample Bag
Material Specification Sheet
B01445WA**

Size: 10" x 15", 25.4 x 38 cm

Capacity: 92 oz, 2,721 ml.

Film: Blend of linear low and low-density virgin polyethylene, 4.0 mils, .102 mm thick, FDA approved for food contact.

Write-On Strip: Horizontal band of thermoset white ink printed on film, compliant with CONEG resolution for heavy metals and European Union Directive for Hazardous Substances.

Closing Tab: One round and one flat aluminum wires covered with plastic tape by a patented process that also covers the wire ends, reduces bag puncture, forms leak proof seal when closed.

Pull Tabs: Made of Tyvek®.

Sterilization: Bag is sterilized with ethylene oxide gas; documentation of sterility is available upon request.

ISO Certification: Manufactured under a quality management system certified to **ISO9001**.

Packaging: 250 bags/box in a liner bag twist tied shut.

Conditions of Use: Not to be used at temperatures above 180F, 82C.