

Date: 10/09/2023

From: Whirl-Pak®

Subject: Whirl-Pak® Sample Bag Material Specification B01365 Stand-Up Bag

Size: 4.5" W x 9" L, 11.5 x 23 cm.

Capacity: 18 oz, 532 ml.

Film: Blend of linear low and low density virgin polyethylene, 3.0 mil, .076 mm thick, FDA approved for food contact.

Construction: A 1", 2.54 cm bottom gusset makes the bag self-standing.

Closing Tab: Two round aluminum wires covered with plastic tape by a patented process that also covers the wire ends, reduces bag puncture, forms leak-proof seal when closed.

Pull Tabs: Made of Tyvek®.

Write-On Strip: Horizontal band of thermoset white ink printed on film, compliant with CONEG resolution for heavy metals and European Union Directive for hazardous substances.

Sterilization: Sterilized with ethylene oxide gas, documentation of sterility is available upon request.

ISO Certification: Manufactured under a quality management system certified to **ISO9001 & ISO13485**.

Packaging: 500 bags/box in a liner bag twist-tied shut.

Condition of Use: Not to be used at temperatures above 180F, 82C.

If you have any questions, or need assistance, please email our Customer Support team at contact@whirl-pak.com. For Frequently asked questions please visit our website www.whirl-pak.com/frequently-asked-questions