



Celebrating over
60 years of service

WORLDWIDE SERVICE TO EDUCATION • HEALTH • AGRICULTURE • INDUSTRY

901 JANESVILLE AVENUE FORT ATKINSON, WISCONSIN 53538-0901
PHONE 920-563-2446 FAX 920-563-8296
www.eNASCO.com E-MAIL info@eNASCO.com

Whirl-Pak® Filter Bag Specification Sheet
B01318WA, B01348WA, B01385WA, B01416WA

Size:	B01318WA 7 ½" W x 12" L (19 x 30 cm)	B01348WA 6" W x 9" L (15 x 23 cm)	B01385WA 3 ¾" W x 7" L (9 ½ x 18 cm)	B01416WA 7 ½" W x 15" L (19 x 38 cm)
Capacity:	55 oz. (1627 ml)	24 oz. (710 ml)	7 oz. (207 ml)	69 oz. (2041 ml)

Film: Blend of low-linear and low-density virgin polyethylene, 4.0 mil (102 mm) thick; B01385 is 3.0 mil (.076 mm) thick, FDA approved for food contact.

Filter: Finely perforated low-density polyethylene, 1 ¼ mil (.032 mm) thick, hole diameter is .013" (.33 mm) 330 microns.

Closing Tab: Aluminum wires covered with plastic tape by a patented process that also covers the wire ends, reduces bag puncture, forms leak-proof seal when closed.

Pull Tabs: Made of Tyvek®.

Sterilization: Sterilized with ethylene oxide gas, documentation of sterility available upon request.

ISO Certification: Manufactured under a quality system certified to **ISO9001**.

Packaging: 250 bags/box in a liner bag twist-tied shut.

Conditions of Use: Not to be used at temperatures above 180° F (82° C) can be frozen to any temperature.