

Date: 10/09/2023

From: Whirl-Pak®

Subject: Whirl-Pak® Sample Bag Material Specification B01451 Stand-Up Bag

Size: 7.5" x 15", 19 x 38 cm.

Capacity: 69 oz, 2,041 ml.

Film: Blend of linear low and low-density virgin polyethylene, 4.0 mils, .102 mm thick, FDA approved for food contact.

Construction: A 1 3/4", 4.45 cm bottom gusset makes the bag self-standing.

Write-On Strip: Horizontal band of thermoset white ink printed on film, compliant with CONEG resolution for heavy metals and European Union Directive for Hazardous Substances.

Closing Tab: Two round aluminum wires covered with plastic tape by a patented process that also covers the wire ends, reduces bag puncture, forms leak proof seal when closed.

Pull Tabs: Made of Tyvek®.

Sterilization: Bag is sterilized with ethylene oxide gas; documentation of sterility is available upon request.

ISO Certification: Manufactured under a quality management system certified to **ISO9001 & ISO13485**.

Packaging: 250 bags/box in a liner bag twist tied shut.

Conditions of Use: Not to be used at temperatures above 180F, 82C.

If you have any questions, or need assistance, please email our Customer Support team at contact@whirl-pak.com. For Frequently asked questions please visit our website www.whirl-pak.com/frequently-asked-questions