

• Diversified Biotech •

Product Specification Sheet

Product Name: Cryo-Tags® and Cryo-Babies® (on rolls)

Face stock: Vinyl – 3.5 mil (100 microns)

Adhesive: Acrylic – 1.0 mil (25 microns)

Liner: 3.2 mil (77.5 microns)

Temperature Resistance: -196°C to 79°C

Adhesion Test: 180° peel strength, 1” wide sample

Application Temperature: Room temperature, 50% relative humidity or less. Surface area must be free of moisture.

Adhesion Time	10 min dwell	3 days at 72°F (22°C)
	at RT - oz. /in.	50% R.H. - oz. /in.
Substrate		
Polypropylene	36	75

Other Tests:

- 1) Submersed in water at 66°C for 100 Hours – no visual change
- 2) Liquid and Vapor Phase Nitrogen – subjected to multiple freeze and thaw cycles – no visual change

Shelf Life: One year from date of manufacture. Store at 72°F and 50% relative humidity.

• Diversified Biotech •

EU REACH

This product is an article, without intended release of a chemical substance, under the Regulation No 1907/2006 of the European Parliament and the Council concerning the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH) (refer to REACH, Article 3(3)). It is not a chemical preparation. Therefore, it is not subject to the pre-registration or the registration process. It does not require a safety data sheet.

This product does not contain at greater than 0.1% by weight a Substance of Very High Concern (SVHC) substance identified according to Article 59 of REACH. This declaration reflects the substances on the candidate SVHC list, effective December 2015

DOMINIQUE DUTSCHER SAS