

TECHNICAL DATA SHEET

Part Number	054818
Description	Microtubes 0.6 ml



Benefits :

The MaxyClear Snaplock Microcentrifuge Tubes all have a frosted cap surface for labeling and frosted panels on the sides for writing.

They also include a piercing port in the center of the cap.

These 99.9%-pure virgin polypropylene microtubes are manufactured using highly polished molds for clarity.

Large writing surface and cap surface for easy labeling; 100 μ L graduations for reference. Maximum Recovery available for critical molecular samples.

Technical specifications :

- ✓ Size - 0.6 mL
- ✓ Non-Sterile
- ✓ Polypropylene
- ✓ Maximum Recovery Surface Treatment
- ✓ Color - Clear

Sales unit :

500 Microtubes
