

# WHEATON® CryoFile® Tissue Storage Box, White

- Made from water repellent fiberboard
- Partitions numbered from 1 – 81 for easy content identification
- Numbering system printed on lid and bottom of box
- Six colors provide easy sample identification
- Water repellent allows for longer durability
- For use with vapor phase of liquid nitrogen
- Directional holes in bottom allow for drainage and orientation of bottom to top of box
- Dimensions (L x W x H): CryoFile® (130 x 130 x 53mm), CryoFile® XL (130 x 130 x 97mm), CryoFile® Tissue Box (130 x 130 x 26mm)
- CryoFile® Tissue Box, holds 25 tissue vials

Catalog Number: **W651610-W**

## Product Overview

Storage box for use with cryogenic vials

## Specifications

Remarks	Type: Fits 5ml tissue vial
Color	White
Colour	White
Quantity	15

*This information is accurate as of 23/05/2024 and is subject to change.*