

## TECHNICAL DATA SHEET

<b>Part Number</b>	055133
<b>Description</b>	Nunc MaxiSorp Plate



### Benefits :

- ✓ Nunc Immunology Plates MaxiSorp Nunc plates
- ✓ 96 wells
- ✓ Flat bottom
- ✓ Well volume 400  $\mu$ l
- ✓ Material : polystyrene
- ✓ Dimension of the plate : 128 mm  $\times$  86 mm
- ✓ Strong affinity for proteins: 600 to 650 ng of IgG per  $\text{cm}^2$
- ✓ Hydrophilic
- ✓ Alphanumeric indications
- ✓ Supplied without lid
- ✓ Non-sterile
- ✓ Delivered with certificate

### Technical characteristics :

- ✓ 96 wells (8x12)
- ✓ Volume per well: 0.4 ml / volume
- ✓ Depth of wells: 11.2 mm

### Sales unit :

60 Plates