



## Thermo Scientific Laboratory Consumables Catalog

Cell Culture | Diagnostics  
Environmental Sample Containers  
Filtration | General Labware | IVF  
Microplates for Drug Discovery | Packaging  
Sample Collection and Processing  
Sample Storage | Vaccines & Biologics

**Quality, Choice**  
**and Innovation**  
for your applications

**Thermo**  
SCIENTIFIC

## ► Thermo Scientific Nunc Centrifuge Tubes, 10 to 11 mL



**Nunc Centrifuge Tubes are tightness tested and are available in polycarbonate and polystyrene.**

### details

- Made in a variety of materials and formats for greater flexibility
- Different cap types and bottom shapes to meet all needs
- Fit most standard centrifuges
- Relative centrifugal force tested (RCF)

**Certifications:** Polystyrene tubes are CE marked.

### Nunc Centrifuge Tubes, 10 to 11 mL

Cat. No.	Description	Shape	Volume	L x W, mm	Max. Force	Sterile	Writing Area	No. per Case
347759	Polycarbonate tube with polypropylene rim-seal cap	Conical	11 mL	110 x 16	6000	Yes	No	480
347708	Polycarbonate tube with polypropylene rim-seal cap	Conical	11 mL	110 x 16	6000	Yes	Yes	480
362308	Polystyrene tube with plug seal polyethylene cap	Flat	10 mL	100 x 16	3000	No	Yes	480
362707	Polystyrene tube with plug seal polyethylene cap	Flat	10 mL	100 x 16	3000	Yes	Yes	480
348224	Polystyrene tube with plug seal polyethylene cap	Round	10 mL	110 x 16	3000	Yes	No	480
347856	Polystyrene tube with plug seal polyethylene cap	Conical	11 mL	110 x 16	3000	Yes	No	480
347880	Polystyrene tube with plug seal polyethylene cap	Conical	11 mL	110 x 16	3000	Yes	Yes	480
337846	Polystyrene tube with plug seal polyethylene cap	Conical	11 mL	110 x 16	3000	Yes	Yes	480