


CATALOG NUMBER	SEE BELOW
WORKING RANGE (µL)	50 - 250µL
MATERIAL	PS
A. BASE LENGTH	MM 127.8 INCH 5.03
B. BASE WIDTH	85.5 3.37
C. TOP LENGTH	123.3 4.84
D. TOP WIDTH	81.2 3.20
E. FLANGE LONG SIDE	7.5 0.3
F. OVERALL HEIGHT	14.4 0.57
G1. DIA TOP	Ø7
G2. DIA BOTTOM	Ø6.2
H. A - I LOCATION (X)	14.3 0.56
J. A - I LOCATION (Y)	11.3 0.45
K. WELL TO WELL SPACING	9 0.35
L. WELL DEPTH	11.4 0.45
M. BOTTOM OF WELL DISTANCE	3.2 0.13
N. FLANGE SHORT SIDE	2.4 0.1
*P. LID LENGTH	127.4 5.02
*Q. LID WIDTH	85.4 3.36
*R. LID HEIGHT	9 3.36
S. STACKED HEIGHT - NO LID	26.9 1.06
*T. STACKED HEIGHT W/LIDS	31.2 1.23
*U. ASSEMBLY HEIGHT W/LID	15.6 0.61



ThermoFisher SCIENTIFIC

Roskilde Site • Kamstrupvej 90 • P.O.Box 280
DK-4000 Roskilde • Denmark
Tel +45 4631 2000 • Fax +45 4631 2042
E-mail: Info-nunc@thermofisher.com

www.nuncbrand.com

**F96 MICROWELL PLATE,
HIGH FLANGE, PS**

SUBJECT: F96 MICROWELL PLATE, HIGH FLANGE, PS

FORM#

APPROVED BY/DATE:

VALID FROM:

CHANGE ORDER NO.:

VALID FROM:

VERSION:

PAGE TYPE: SCALE:

MATERIAL:

DRAWING NO.:

ITEM NUMBER:

VERSION:

SCALE:

META INFORMATION:

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC

CATALOG NO	SURFACE	COLOR	STERILE	WITH LID*	BAR CODE	UNITS/ PK/CASE	CATALOG NO	SURFACE	COLOR	STERILE	WITH LID*	BAR CODE	UNITS/ PK/CASE
156545	CELL CULTURE	CLEAR	+	-	-	10/180	260844	NON-TREATED	CLEAR	-	-	+	10/180
161093	CELL CULTURE	CLEAR	+	-	-	10/180	260860	NON-TREATED	CLEAR	+	+	-	10/180
167008	CELL CULTURE	CLEAR	+	-	-	1/50	260887	NON-TREATED	CLEAR	+	+	+	10/180
168055	CELL CULTURE	CLEAR	+	-	-	1/50	260895	NON-TREATED	CLEAR	-	+	-	10/180
152038	COLLAGEN	CLEAR	-	-	-	5/20	456529	POLYSORP	CLEAR	-	-	-	10/180
152039	POLY-D- LYSINE	CLEAR	-	-	-	5/20	456537	MAXISORP	CLEAR	-	-	-	10/180
176747	CC™	CLEAR	+	-	-	5/50	460984	MAXISORP	CLEAR	-	-	+	10/180
256510	NON-TREATED	CLEAR	+	-	-	10/180	467320	MEDISORP	CLEAR	-	-	-	5/60
260636	NON-TREATED	CLEAR	-	-	-	10/180	467340	MULTISORP	CLEAR	-	-	-	5/60

* DIM'S APPLY ONLY WHERE "WITH LID" = +. USING NUNC LID 2205

NOTES:

- ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
- THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATION PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.