

TECHNICAL DATA SHEET

Part number	055730
Description	Ramboldi® 30ml Polystyrene Universals



Description :

- ✓ In polystyrene
- ✓ Cap material : polypropylene
- ✓ Fitted with proof, flow seal cap
- ✓ Supplied sterile by EtO
- ✓ Capacity : 30 ml
- ✓ Label : plain
- ✓ Dimensions : 30.5 x 93.5 mm
- ✓ Inner packs quantity : 50

Pieces by sales :

400 Tubes
