

WHEATON® CENTRIFUGE TUBES

Ideal for standard laboratory sample management
and centrifugation needs.



DWK LIFE SCIENCES
DURAN WHEATON KIMBLE

WHEATON® POLYPROPYLENE TUBES ARE IDEAL FOR STANDARD LABORATORY CENTRIFUGATION NEEDS.

Manufactured from ultra-clear virgin polypropylene, the WHEATON® Centrifuge Tube range meets the needs of most laboratory users for routine centrifugation applications. Available in sizes 15ml and 50ml, with a leak-free cap and large writing area, they are ideal for a wide range of sample handling applications.

CENTRIFUGE TUBE FEATURES

- Manufactured from virgin ultra-clear polypropylene for good sample visibility
- Flat cap design:
 - Offers excellent leak-free performance
 - Easy to write on and mark samples
- Solvent-resistant printed graduations with large white marking strip for labelling
- RCF values up to 9500xg
- Non-pyrogenic to < 02. 5EU/mL
- RNase and DNase free
- Packaging options:
 - Packed in a polystyrene rack
 - Bulk packed in a peel open-bag
- Supplied sterile
- Operating temperatures -140°C to 121°C
- Easy to handle with gloved hands



PEEL BAG PACK



FOAM RACK PACK



ORDERING INFORMATION

| Code | Description | Max RCF (xg) | Graduation Range (ml) | Graduation Intervals (ml) | Sterility | Inner Case Qty. | Outer Case Qty. |
|--------|---|--------------|-----------------------|---------------------------|-----------|-----------------|-----------------|
| 380220 | 15ml Centrifuge Tube, Polypropylene, Bagged, Sterile | 9000 | 2 to 14.5 | 0.5 | ETO | 50 | 500 |
| 380230 | 15ml Centrifuge Tube, Polypropylene, Racked, Sterile | 9000 | 2 to 14.5 | 0.5 | ETO | 15 | 500 |
| 380100 | 50ml Centrifuge Tube, Polypropylene, Bagged, Sterile | 9500 | 5 to 47.5 | 2.5 | ETO | 25 | 500 |
| 380110 | 50ml Centrifuge Tube, Polypropylene, Racked, Sterile | 9500 | 5 to 47.5 | 2.5 | ETO | 25 | 500 |
| 380240 | 50ml Centrifuge Tube, Polypropylene, Self-Standing, Bagged, Sterile | 7000 | 5 to 47.5 | 2.5 | ETO | 25 | 500 |

WH.CT-ENG-V7-29.07.2020



DWK Life Sciences Limited
Unit 4, Riverside 2, Campbell Road,
Stoke-on-Trent, ST4 4RJ,
UK

Phone: +44 (0)1782 444406
Fax: +44 (0)1782 940436
enquiries.uk@dwk.com
www.dwkltd.com