


# WHEATON®

## PLASTIC PASTEUR PIPETTES

Ideal for the accurate volume transfer of liquid and bacteriological samples.



**DWK LIFE SCIENCES**  
DURAN WHEATON KIMBLE

# WHEATON® SINGLE-USE POLYETHYLENE PASTEUR PIPETTES ARE IDEAL FOR GENERAL LIQUID TRANSFER APPLICATIONS AND EXTRACTING SAMPLES FOR ANALYSIS.

WHEATON® polyethylene pasteur pipettes are an everyday laboratory essential. Available in a range of sizes from 1ml to jumbo-sized pipettes with fine tip variants, there is a product to suit most user applications for use in areas such as blood banking, bacteriology, general life science and industry. The integral bulb design and single-use feature eliminate the need for plugging or autoclaving, helping to reduce the potential of cross-contamination during everyday sample handling.

## PASTEUR PIPETTE FEATURES

- Manufactured from soft, flexible, inert non-toxic polyethylene
- Ideal for safe and easy transfer of liquid samples
- Clear graduations on the 1ml and 3ml versions support ease of volume control
- Suitable for use as a freeze vial when the bulb is fully heat sealed
- Jumbo pipettes are ideal for larger volume transfer or accessing liquids in taller vessels
- Easy to handle with gloved hands
- Non-toxic material can be safely incinerated for disposal
- Available as both sterile and non-sterile options in easy-to-open, paper-plastic peel pouches



## ORDERING INFORMATION

Code	Description	Sterility	Inner Case Qty.	Outer Case Qty.
711115	1ml Graduated Pasteur Pipette, Bulk	NS	500	3000
726129	1ml Graduated Pasteur Pipette, Individual Wrap	IRR	1	500
711117S/10	1ml Graduated Pasteur Pipette, Pack of 10	IRR	10	1500
711117S/20	1ml Graduated Pasteur Pipette, Pack of 20	IRR	20	1500
711116	3ml Graduated Pasteur Pipette, Bulk	NS	500	3000
726128	3ml Graduated Pasteur Pipette, Individual Wrap	IRR	1	500
711118S/10	3ml Graduated Pasteur Pipette, Pack of 10	IRR	10	1500
711118S/20	3ml Graduated Pasteur Pipette, Pack of 20	IRR	20	1500
711117F	1ml Fine-Tip Pasteur Pipette, Bulk	NS	500	3000
711117F1S	1ml Fine-Tip Pasteur Pipette, Individual Wrap	IRR	1	500
711117F10S	1ml Fine-Tip Pasteur Pipette, Pack of 10	IRR	10	1500
711117F20S	1ml Fine-Tip Pasteur Pipette, Pack of 20	IRR	20	1500
711230	Jumbo Pasteur Pipette	NS	100	600
711117MF/3000	Fine-Tip Micro-Pasteur Pipette, Bulk	NS	500	3000

WH.PPIP-ENG-V5-02.07.2020



DWK Life Sciences Limited  
Unit 4, Riverside 2, Campbell Road,  
Stoke-on-Trent, ST4 4RJ,  
UK

Phone: +44 (0)1782 444406  
Fax: +44 (0)1782 940436  
enquiries.uk@dwk.com  
www.dwkltd.com