



WHEATON® 30ML POLYPROPYLENE UNIVERSALS

Versatile container range

Ideal for your sample collection, transport and storage needs

- Ideal for general patient sample collection, storage and transport needs
- Manufactured from virgin, medical grade, polystyrene
- Fitted with leakproof flow seal cap
- Available plain, pre-filled with boric acid for the preservation of urine samples or with an integral spoon for fecal sample collection
- Supplied sterile by EtO
- Convenient inner bags of 50 pieces within each case
- CE Marked in accordance with European Directive 98/79/EC



Product Specification

Item Code	Label	Inner Pack Qty	Cap Material	Description	Sterile	Diameter x height mm	Capacity ml	Pack quantity
327150P P	Plain	50	Polypropylene		ETO	30.5 x 93.5	30	400
327151P P	No	50	Polypropylene	30ml Polypropylene Universal, No Label, Sterile	ETO	30.5 x 93.5	30	400
327152P P	Printed	50	Polypropylene	30ml Polypropylene Universal, Printed Label, Sterile	ETO	30.5 x 93.5	30	400
327153P P	Printed	50	Polypropylene	30ml Polypropylene Universal, Boric Acid, Sterile	ETO	30.5 x 93.5	30	400
327155P P	Printed	50	Polypropylene	30ml Polypropylene Universal, Spoon, Sterile	ETO	30.5 x 93.5	30	400

DOMINIQUE DUTSCHER SAS