




CATALOG NUMBER.....	159933	159934		
RECOMMENDED WORKING VOL.:	150	150		
MAXIMUM CAPACITY IN UPRIGHT POSITION (ml).....	800	800		
CAP TYPE.....	VENT/CLOSE	FILTER		
CULTURE AREA (cm ²).....	225	225		
MATERIAL	PS	PS		
FLASK	PS	PS		
CAP	PS	PS		
	BLUE HDPE	BLUE HDPE		
	INCH	MM	INCH	MM
A (NECK I.D.).....	1.04	26.4	1.04	26.4
B (LENGTH W/O CAP).....	9.06	230.1	9.06	230.1
C (LENGTH W/ CAP).....	9.13	231.9	9.13	231.9
D (WIDTH).....	5.54	140.7	5.54	140.7
E (HEIGHT).....	1.57	39.9	1.57	39.9


1. ALL DIMENSIONS ARE NOMINAL AND IN INCHES AND MILLIMETERS UNLESS OTHERWISE NOTED.
2. THIS IS NOT A PRODUCTION DRAWING, BUT IS FOR INFORMATIONAL PURPOSES ONLY, AND IS SUBJECT TO CHANGE WITHOUT NOTICE.



SIDE VIEW



TOP VIEW



END VIEW OF 159933



END VIEW OF 159934

SPECS ON BACK

NOT TO SCALE

CONFIDENTIAL

THIS PRODUCT SPECIFICATION SHEET CONTAINS CONFIDENTIAL AND/OR PROPRIETARY INFORMATION AND IS TO BE USED ONLY FOR THE EXPRESSED PURPOSES AGREED UPON BY AN AUTHORIZED OFFICER OF THERMO FISHER SCIENTIFIC.

ThermoFisher
SCIENTIFIC

Thermo Fisher Scientific
75 Panorama Creek Drive
Rochester, NY 14625-2385, USA

NAME	MC ²	MP ²	PVC	REV.	0306
DATE	MP	MP	3/9/06		
DESCRIPTION:	PSS 15993X				

EASY FLASK 225CM² WITH CLOSURE