



Thermo Scientific
Nunc Serological Pipettes

research with confidence

Thermo
SCIENTIFIC

Thermo Scientific Nunc Serological Pipettes

Our newly expanded portfolio of serological pipettes features a full range of size and packaging options—offering the quality, accuracy and performance today's scientists require during every stage of cell culture research.

We are dedicated to meeting the unique demands of your research. At every stage of cell culture, our new selection of serological pipettes offers the quality, accuracy, size versatility and packaging options you need to advance your research with complete confidence.

Cell Isolation
and Preparation

Characterization
and Analysis

Harvest

Storage

Growth
and Passage

Serological Portfolio Highlights

- Use new, sleek, slotted chelates for easy measuring under the hood
- Choose 10 pipette sizes ranging from 1 mL to 50 mL with nose-clipper tips and fully finished graduations for precision
- High-precision chelates for higher volume use with less packaging waste
- Delicate iron chelate wrapped chelates in either paper/kraft peel or plastic film
- Free fluoride ions with new chelates
- Clear, 100% virgin polycarbonate

- Accurate, consistent pipette pipettes plugged and sealed
- Delicate iron chelate wrapped chelates
- 100% virgin polycarbonate to maximize clarity
- Coated tips to optimize flow
- Delicate iron chelates with double-walled soft mouthpiece
- Nose-clipper tips for stable precision
- Coated packaging for ease in sorting and selecting the correct size

- Packaging options that accommodate every lab
- Bulk packs, paper/kraft peel or plastic film and in pre-filled or for easy opening
- Bulk packs for frequent use reduce packaging waste
- Multiple sizes for use in a wide range of applications



Thermo Scientific Nunc Serological Pipettes

Ordering Information

Cat No	Description	Unit per bag	Unit per case
Individually wrapped, paper peel packaging, plugged			
170353	1 ml pipette	100	1000
170354	2 ml pipette	100	500
170355	5 ml pipette	60	200
170356	10 ml pipette	50	200
170357	25 ml pipette	50	200
170358	50 ml pipette	25	100
Shortie, Individually wrapped paper/plastic peel packaging, plugged			
170360	5 ml shortie pipette	50	200
170361	10 ml shortie pipette	50	200
Individually wrapped, plastic film packaging, plugged			
170364	1 ml pipette	100	1000
170365	2 ml pipette	100	500
170366	5 ml pipette	50	200
170367	10 ml pipette	50	200
170368	25 ml pipette	50	200
170369	50 ml pipette	25	100
Pipettes, Bulk Packed, Plugged			
170371	1 ml pipette	50	1000
170372	2 ml pipette	50	500
170373	5 ml pipette	25	500
170374	10 ml pipette	25	500
170375	25 ml pipette	25	200
170376	50 ml pipette	25	100



Thermo Scientific S1 Pipet Filler

A lightweight, cordless pipet filler that allows longer, fatigue-free pipetting—an excellent companion to our serological pipettes.

For more information, go to
www.thermoscientific.com/s1pipetfiller

© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trade-marks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

www.thermoscientific.com

ANZ: Australia: 1300 725 232, New Zealand: 0800 933 066

Asia: China Tel: +86 21 50451118 or 4004504511; India Tel: +91 22 5716 2210; India Toll-Free: 1 800 22 8934;

Japan: 011 31006 1611; Other Asian countries: 16 60271212

Europe: Austria: +43 1 801 40 0; Belgium: +32 51 73 42 41; Denmark: +45 4631 2000; France: +33 2 2803 21 80;

Germany: +49 6 34 90 89 40; Germany Toll-Free: 0800 1 536 316; Italy: +39 02 95039 552 or 434 254 475;

Netherlands: +31 76 521 444; Nordic/Baltic countries: +358 9 329 100; Russia/CIS: +7 (812) 703 02 15;

Spain/Portugal: +34 93 223 09 18; Switzerland: +41 44 454 12 12; UK/Ireland: +44 870 609 8203

North America: USA/Canada: 1 800 541 0000, USA Toll-Free: 800 671 4527

South America: USA sales support: +1 985 869 7198

Countries not listed: +49 6194 90 8910 or +33 2 2803 2182

BLS/SEPP/0412

Thermo
SCIENTIFIC

Part of Thermo Fisher Scientific