



Thermo Scientific
Nunc Serological Pipettes

research with confidence

DOMINIQUE DUTOUR SAC

Thermo
SCIENTIFIC

Thermo Scientific Nunc Serological Pipettes

Our newly expanded portfolio of serological pipettes features a full range of size and packaging options – offering the quality, accuracy and performance today's scientists require during every stage of cell culture research.

Cell isolation and preparation



Characterization and Analysis



Harvest



Storage



We are dedicated to meeting the rigorous demands of your research, at every stage of cell culture. Our new collection of serological pipettes offers the quality assurance, size variability and packaging options you need to advance your research with complete confidence.

Growth and passage



Serological Portfolio Highlights


- Use new sleek shortie pipettes for easy maneuvering under the hood
- Choose multiple sizes, ranging from 1 mL to 50 mL, all with non-dripping tips and clearly marked graduations for precision
- Pick bulk packed options for higher volume use with less packaging waste
- Select individually wrapped pipettes in either paper/plastic peel or plastic film
- See the difference with new ultra clear, 100% virgin polystyrene

Accurate, disposable plastic pipettes plugged and sterilized

- Sterility assurance level (SAL) of 10^{-6}
- 100% virgin polystyrene for maximum clarity
- Certified non-pyrogenic
- Bold clear graduations with easy-to-read scale markings
- Non-dripping tip for pipette precision
- Color-coded packaging for ease in sorting and selecting the correct size

Packaging options that accommodate every lab!

- Individually wrapped; paper/plastic or plastic film; peel or pop through for easy opening
- Bulk packed for frequent use; reduces packaging waste
- Multiple sizes for use in a wide range of applications



DOMINIQUE DUTCH

Thermo Scientific Nunc Serological Pipettes

Ordering Information

| Cat No | Description | Unit per bag | Unit per case |
|---|-----------------------|--------------|---------------|
| Individually wrapped, paper peel packaging, plugged | | | |
| 170353 | 1 ml pipette | 100 | 1000 |
| 170354 | 2 ml pipette | 100 | 500 |
| 170355 | 5 ml pipette | 50 | 200 |
| 170356 | 10 ml pipette | 50 | 200 |
| 170357 | 25 ml pipette | 50 | 200 |
| 170358 | 50 ml pipette | 25 | 100 |
| Shortie, Individually wrapped paper/plastic peel packaging, plugged | | | |
| 170360 | 5 ml shortie pipette | 50 | 200 |
| 170361 | 10 ml shortie pipette | 50 | 200 |
| Individually wrapped, plastic film packaging, plugged | | | |
| 170364 | 1 ml pipette | 100 | 1000 |
| 170365 | 2 ml pipette | 100 | 500 |
| 170366 | 5 ml pipette | 50 | 200 |
| 170367 | 10 ml pipette | 50 | 200 |
| 170368 | 25 ml pipette | 50 | 200 |
| 170369 | 50 ml pipette | 25 | 100 |
| Pipettes, Bulk Packed, Plugged | | | |
| 170371 | 1 ml pipette | 50 | 1000 |
| 170372 | 2 ml pipette | 50 | 500 |
| 170373 | 5 ml pipette | 25 | 500 |
| 170374 | 10 ml pipette | 25 | 500 |
| 170375 | 25 ml pipette | 25 | 200 |
| 170376 | 50 ml pipette | 25 | 100 |



Thermo Scientific S1 Pipet Filler

A lightweight, cordless pipet filler that allows longer, fatigue-free pipetting—an excellent companion to our serological pipettes.

For more information, go to www.thermoscientific.com/s1pipetfiller

© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

www.thermoscientific.com

ANZ: Australia: 1300 735 292, New Zealand: 0800 933 966

Asia: China Toll-free: 800-810-5118 or 400-650-5118; India: +91 22 6716 2200, India Toll-free: 1 800 22 8374; Japan: +81 3 5826 1616; Other Asian countries: 65 68729717

Europe: Austria: +43 1 801 40 0; Belgium: +32 53 73 42 41; Denmark: +45 4631 2000; France: +33 2 2803 2180; Germany: +49 6184 90 6940, Germany Toll-free: 08001-536 376; Italy: +39 02 95059 552 or 434-254-375; Netherlands: +31 76 571 4440; Nordic/Baltic countries: +358 9 329 100; Russia/CIS: +7 (812) 703 42 15; Spain/Portugal: +34 93 223 09 18; Switzerland: +41 44 454 12 12; UK/Ireland: +44 870 609 9203

North America: USA/Canada +1 585 586 8800; USA Toll-free: 800 625 4327

South America: USA sales support: +1 585 899 7198

Countries not listed: +49 6184 90 6940 or +33 2 2803 2180

BRLSPSERPIP 0412

Thermo
SCIENTIFIC

Part of Thermo Fisher Scientific