



# PRODUCT DATA SHEET

**COWIE TECHNOLOGY GROUP LTD**  
RIDGEWAY, COULBY NEWHAM  
MIDDLESBROUGH, CLEVELAND  
TS8 0TQ, ENGLAND

TEL: +44 (0)1642 599190  
FAX: +44 (0)1642 596810

enquiries@cowie.com  
www.cowie.com

Item: Homogeniser - Glass Vessel	
Product Code/Ref: 011.115	
No. Pages:	1

Material: Borosilicate glass.

Capacity: 15ml.

Bore: 19mm.

Height: 155

Sterilise: Thermal, chemical gamma radiation.

GRC 01/16

PTFE used in the manufacture of this product is VIRGIN PTFE GRADE TF1620 manufactured by 3M® DYNEON®

**3M® DYNEON® TF1620 – COMPLIANCE STATUS**

- FDA – Directive FDA 21 CFR 177.1550
- USP Class VI – USP 27, NF 22, 2004 for Plastics at 70°C.
- RoHS 2 – Directive 2011/65/EU & Dangerous Substances – 76/769/EC
- EU REACH Regulations – SVHC Current Candidate List
- Bisphenol A – Absent
- Latex – Absent
- ADCF – Animal Derived Components Free (BSE/TSE). GMO Free
- Endotoxins – Free
- POP, PFOS, PFOA, APFO, DEHIP, VOC Free
- ODC – Class 1
- Dodd-Frank Act – Gold, Coltan, Cassiterite, Wolframite and their derivatives originating from the Democratic Republic of Congo - Absent

It is the responsibility of the user to ensure that the products are suitable for their intended use

