



# PRODUCT DATA SHEET

**COWIE TECHNOLOGY GROUP LTD**  
**RIDGEWAY, COULBY NEWHAM**  
**MIDDLESBROUGH, CLEVELAND**  
**TS8 0TQ, ENGLAND**

**TEL: +44 (0)1642 599190**  
**FAX: +44 (0)1642 596810**

**enquiries@cowie.com**  
**www.cowie.com**

Item: Homogeniser - Plunger	
Product Code/Ref: 011.315	
No. Pages:	1

Material: PTFE / Type 304 S/S Shaft  
Capacity: 15ml.  
Bore: 19mm.  
Height: 270mm.  
Sterilise: Thermal, chemical NOT gamma radiation.

GRC 01/16

PTFE used in the manufacture of this product is VIRGIN PTFE GRADE TF1620 manufactured by 3M® DYNEON®

**3M® DYNEON® TF1620 – COMPLIANCE STATUS**

- FDA – Directive FDA 21 CFR 177.1550
- USP Class VI – USP 27, NF 22, 2004 for Plastics at 70°C.
- RoHS 2 – Directive 2011/65/EU & Dangerous Substances – 76/769/EC
- EU REACH Regulations – SVHC Current Candidate List
- Bisphenol A – Absent
- Latex – Absent
- ADCF – Animal Derived Components Free (BSE/TSE). GMO Free
- Endotoxins – Free
- POP, PFOS, PFOA, APFO, DEHIP, VOC Free
- ODC – Class 1
- Dodd-Frank Act – Gold, Coltan, Cassiterite, Wolframite and their derivatives originating from the Democratic Republic of Congo - Absent

It is the responsibility of the user to ensure that the products are suitable for their intended use

