

TECHNICAL DATA SHEET

Part Number	059918
Description	Aniosyme X3 for pre-disinfection of instrumentation



Benefits :

- ✓ Reinforced cleaning and pre-disinfection of medico-surgical instrumentation, medical devices and endoscopic equipment
- ✓ Very high performing detergency
- ✓ Cleaning and degreasing action
- ✓ Clear blue liquid
- ✓ optimised antimicrobial efficacy from 5 minutes
- ✓ Cleaning in ultrasonic bins
- ✓ Collection of contaminated instrumentation
- ✓ Tri-enzyme complex combined with a new generation onionic surfactant
- ✓ Neutral pH : compatible with aluminium alloys
- ✓ Chloride free formula : non corrosive on materials
- ✓ pH of pure product : around 8
- ✓ 200 doses, 25 ml each

Sales unit :

200 doses
