

## 2- sliding door tall safety cabinet - 2 compartments

### COMPLIANCE

### ADVANTAGES

- Made in France
- INRS recommendations compliant
- Storage of flammables, corrosives and toxics

### BUT ALSO...

- Compliant to the Public Health Code articles R. 5132-66 and R. 5132-68
- INRS recommendations compliant
- Holes for the connection to the ventilation
- Normalized labelling provided with the equipment
- Allows the separation of different types of chemicals
- Locking system



### TECHNICAL FEATURES

Model	<b>A2T</b>	Shipping dimensions HxWxD (mm)	<b>2025 x 1200 x 800</b>
Cabinet materials	<b>10/10th epoxy painted steel structure</b>	Ventilation	<b>Inlet located at the bottom of the cabinet, outlet located at the top</b>
Doors color	<b>white RAL 9010</b>	Openings diameter (mm)	<b>100</b>
Number of doors	<b>2</b>	Number of compartments	<b>2</b>
Doors type	<b>swing doors</b>	HS Code	<b>94032080</b>
Door closure	<b>manual closing</b>		

### ACCESSOIRES COMPLÉMENTAIRES CONSEILLÉS - CAISSONS DE FILTRATION OU DE VENTILATION

- Filtration box with activated carbon filter for organic and corrosive vapors (Réf : CDF-ACORG)
- Filtration box built in PVC with activated carbon filter for organic and corrosive vapors (Réf : CDF-P)
- Ventilation box for a connection to the outside (Réf : CDV-A)
- Recycling air box Labopur - NF X 15-211 (Réf : H50C)
- Connecting kit (Réf : KRC)

### OPTIONS

- Retention shelf for A1T - A2T - A3T - A4T (Réf : ES1T)

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Storage capacity (L)	Containment sump	Shelves	Load capacity	Weight (kg)
A2T	1825 x 1120 x 545	2x (1800 x 550 x 475)	240		8	50	160.00