

MF25



Wall-mounted medicine cabinet - 4 shelves

COMPLIANCE

ADVANTAGES

- Locking system
- Storage of pharmaceuticals products
- Storage of first-aid equipments



TECHNICAL FEATURES

Model	MF25	Hinges	right
Cabinet materials	10/10th epoxy painted steel structure	Door closure	manual closing
Body color	white RAL 9010	Shipping dimensions HxWxD (mm)	780 x 600 x 800
Doors color	white RAL 9010	Number of compartments	1
Number of doors	1	HS Code	94032080
Doors type	swing doors		

Model	External dimensions HxWxD (mm)	Internal dimensions HxWxD (mm)	Eye-washer flow (L/min)	Shower flow (L/min)	Weight (kg)
MF25	580 x 470 x 205	570 x 460 x 170	0	0	15.00



TRIONYX
1^{er} FABRICANT FRANÇAIS
Armoires de sécurité pour produits dangereux

138 avenue Henri BARBUSSE 92140 Clamart (France)
Tél. 01 46 45 80 00 - Fax : 01 46 45 85 98
www.trionyx.fr - trionyx[at]trionyx.fr