

ORG50

NF X 15-211  

Filtro ha carbone attivo per vapori organici - for filtration box H50C (or old cabinets range 12.X before 2002)

CONFORMIDAD

ADVANTAGES

- Specific filtration of the chemicals
- NF X 15-211 compliant
- Storage of flammables, corrosives and toxics



CARACTERISTICAS TECNICAS

Modelo	ORG50			Código Aduanero	38021000		
Modelo	Dimensiones exteriores alt x anc x pro (mm)	Dimensiones interiores alt x anc x pro (mm)	Capacidad de almacenamiento (L)	Cubetas de retención	Estanterías	Carga admisible	Peso (kg)
ORG50	70 x 390 x 390		0		0	0	8.00