

# BioMixer™ and MiniMixer™ Benchmark 3-Dimensional Shakers



**BioMixer™**  
12 x 12 in. platform

**MiniMixer™**  
10.5 x 7.5 in. platform



B3D1320  
shown with stacking platform

- Nutating motion for thorough, uniform mixing
- Dimpled mat, perfect for blood and centrifuge tubes
- Stacking platforms and flat mats available
- Incubator & cold room safe

## Technical Data

Speed:	24 rpm*
Tilt Angle:	+/- 20 degrees
Platform Size:	
MiniMixer:	10.5x7.5 in. / 27x19 cm
BioMixer:	12x12 in. / 30x30 cm
Load Capacity:	1.6 kg / 3.5 lbs.
Operating Temp. Range:	+4°C to +65°C
Dimensions (WxDxH):	
MiniMixer:	10.5 x 7.5 x 8 in.
BioMixer:	12 x 12 x 8.3 in.
Weight:	2.0 kg / 4.4 lb
Warranty:	2 Years
Electrical:	115V AC, 60 Hz, 0.4A 230V AC, 50 Hz, 0.2A

\*Speed of 230V model is 20 rpm

These platform shakers provide a 3-dimensional “nutating” motion for moderately aggressive mixing. Two platform sizes and two rubber mat configurations allow for thorough, foam-free mixing in blood tubes and centrifuge tubes, as well as in a variety of flat vessels.

Both models come standard with autoclavable dimpled mats for tubes. Optional flat mats are available for mixing other vessel types. Stacking platforms are also available, effectively doubling the usable surface area of each model. Even the smallest single platform configuration has 24% more platform work surface than traditional nutating mixers.

The **MiniMixer™** and **BioMixer™** can handle loads up to 2kg. Both have maintenance free brushless motors and are safe for use in cold rooms and incubators.

## Ordering Information

<b>B3D1020*</b>	<b>MiniMixer, 10.5 x 7.5in. platform with dimpled mat</b>
B3D-STACK-D	Stacking Platform, dimpled mat, 10.5 x 7.5 in.
B3D-FLAT	Flat Mat, 10.5 x 7.5 in.
B3D-DIMPLED	Dimpled Mat, 10.5 x 7.5 in.
<b>B3D1320*</b>	<b>BioMixer, 12 x 12in. platform with dimpled mat</b>
BR1000-STACK-D	Stacking Platform, dimpled mat, 12 x 12 in.
BR1000-FLAT	Flat Mat, 12 x 12 in.
BR1000-DIMPLED	Dimpled Mat, 12 x 12 in.

\* To order a product in 230V, please add -230 to item number.

**Benchmark**   
Scientific

**Benchmark Scientific Inc.**  
PO Box 709  
Edison, NJ 08817  
Phone: 908-769-5555  
Fax: 908-222-1864

Web: [www.BenchmarkScientific.com](http://www.BenchmarkScientific.com)  
Email: [Info@BenchmarkScientific.com](mailto:Info@BenchmarkScientific.com)