

myBLOCK™

Personal sized drybaths



- + Compact, fits on any crowded benchtop
- + Available in heat only or heat & cool versions
- + Exchangeable blocks, for tubes 0.2 to 50ml
- + Simple touchpad, digital control
- + Link 2 consecutive temperatures (cooling model only)

Technical Data

| | BSH200 | BSH300 |
|------------------------|----------------------|----------------------|
| Temp. Range: | RT +5° to 100°C | RT -25° to 100°C |
| Min. Programmable Temp | 20°C | 0°C |
| Temp. Accuracy: | +/- 0.5°C | +/- 0.5°C |
| Temp. Increments: | 0.1°C | 0.1°C |
| Temp. Uniformity: | +/- 0.2°C | +/- 0.2°C |
| Dimensions: | 4.5 x 5.9 x 4 in. | 4.3 x 6.4 x 5.5 in. |
| (w x d x h) | 11.2 x 15 x 10 cm | 11 x 16.2 x 14 cm |
| Weight: | 2 lbs / 1 kg | 2 lbs / 1 kg |
| Electrical: | 100 - 240V, 50-60Hz* | 100 - 240V, 50-60Hz* |

* 12V vehicle power adapter also available

With its compact construction, broad temperature range, and multiple block options, the *myBlock Mini* is the perfect tool for laboratory incubations. The simple touch pad controls and a digital display provide for "set and walk away" temperature selection and unrivaled accuracy.

While perfect for bench use, the optional 12V adapter expands the utility of the *myBlock Mini* to field use. At less than 4.3 inches wide, *myBlock Mini* is truly the first personal block incubator.

Ordering Information:

| | |
|---------------|--|
| BSH200 | MyBlock Mini Dry Bath with block lifter (Blocks sold separately) |
| BSH300 | MyBlock Mini Dry Bath with cooling with block lifter (Blocks sold separately) |
| BSH100-02 | Block, 40 x 0.2ml or 5 x PCR strips |
| BSH100-05 | Block, 24 x 0.5ml |
| BSH100-15 | Block, 15 x 1.5ml |
| BSH100-1520 | Block, 15 x 1.5 or 2.0ml |
| BSH100-150 | Block, 4 x 15ml* |
| BSH100-500 | Block, 2 x 50ml* |
| BSH100-5MT | Block, 6 x 5ml centrifuge tubes* |
| BSH100-1232 | Block, 15 x 2ml autosampler vials |
| BSH100-1213 | Block, 15 x 12/13mm tubes/vacutainers* |
| BSH100-CV | Block, 8 x cuvettes* |
| BSH100-A12 | 12V vehicle power adapter |

*Not compatible with clear cover



Also Available...



myBLOCK™ HL
With Heated Lid

- Prevents evaporation and condensation
- Ideal for concentration sensitive reactions

Item: BSH200-HL*
MyBlock HL with Heated Lid

* For use with 0.2ml, 0.5ml and 1.5ml blocks only.

Benchmark 
Scientific

PO Box 709 - Edison, NJ 08817
Phone: 908-769-5555 - Fax: 732-313-7007
Web: www.BenchmarkScientific.com
Email: Info@BenchmarkScientific.com