

# myFUGE™ 12

## Mini Centrifuge

- + COMBI-Rotor™ for tubes and strips
- + Twice the capacity of traditional mini centrifuges
- + Near silent operation
- + Starts & stops with opening/closing of the lid



### Technical Data

Speed:	5,500rpm / 2,000 xg
Capacity:	12 x 1.5 / 2.0ml tubes 32 x 0.2ml PCR tubes 4 x PCR strips (8x0.2ml)
Dimensions:	5.5 x 7.9 x 4.4 in 14 x 20 x 11.2 cm
Weight:	2.1 lbs / 5 kg
Electrical:	100 to 240V, 50-60Hz

Unlike traditional mini centrifuges, the new MyFuge 12 eliminates the need to change rotors when switching between microtubes and PCR strips. The included, unique COMBI-Rotor is all that is required for running 12 microtubes and 4 PCR strips simultaneously.

With a fixed speed that produces 2,000 xg, this centrifuge is perfect for quick spin downs. Simply close the lid and the unit quickly ramps up to 5500 rpm. Open the lid, and the rotor quickly decelerates for removal of samples.

At just 5.5 inches wide and less than 4.5 inches tall, the MyFuge 12 truly is a personal centrifuge with unmatched capacity and flexibility.



Unique COMBI-Rotor



Compact, low profile design

### Ordering Information

C1012	MyFUGE Mini Centrifuge, with COMBI-Rotor for microtubes and PCR strips/tubes
C1008-A5-6	Adapters, 0.5ml, pack of 6
C1008-A2-6	Adapters, 0.2ml, pack of 6

\*To order in 230V, add (-E) to the item number.

**Benchmark**  
Scientific

**Benchmark Scientific Inc.**  
PO Box 709  
Edison, NJ 08817  
Phone: 908-769-5555  
Fax: 908-222-1864

Web: [www.BenchmarkScientific.com](http://www.BenchmarkScientific.com)  
Email: [Info@BenchmarkScientific.com](mailto:Info@BenchmarkScientific.com)