

# Snap-Cap 5mL MacroTubes®

- Compatible with standard equipment
- Withstands 25,000xg\*
- Easy-open flip top
- Unsurpassed clarity
- Dual graduations

5mL MacroTubes™ take all of the conveniences of standard 1.5mL tubes and scales them up for larger volumes. These polypropylene tubes come with an attached cap that not only provides a tight, leak proof seal, it can also be effortlessly opened (or closed) with one hand. The cap is designed to provide a secure seal from -86°C to +80°C. CapLock Clips™ are available to maintain the integrity of the seal at extreme temperatures. In incubated/heated applications, they may be "boiled" for up to 10 minutes. At temperatures above 80°C, the use of the CapLock Clip™ is recommended.

The tubes are made of maximum clarity, low-binding virgin polypropylene and are certified RNase & DNase free. Sterile tubes are also pyrogen free. They are sized to be compatible with all Benchmark™ (USA) and standard German-designed equipment, including centrifuges, dry bath block heaters, tube racks and thermal shakers/vortexers. The 16mm tube diameter is identical to the diameter of standard 15mL tubes, which ensures further compatibility with centrifuge rotors, etc. In many applications, these tubes have replaced 15mL tubes, dramatically reducing plastic waste, storage space and shipping costs. The potential contamination of pipettor shafts is greatly reduced because the 56mm total tube length (vs. 120mm for standard 15mL tubes) is short enough for standard 1mL or 5mL tips to reach the tube's conical bottom.

\* When properly supported in appropriate rotor.



NEW



**C2590**  
Two-tiered Pipetting Rack  
for 16 x 5mL tubes



**C2586**  
CapLock™ Clip recommended  
for extreme temperatures



**C2592**  
High volume gripper rack  
for 50 x 5mL tubes



**C2581**  
Polycarbonate cryobox for  
25 x 5mL tubes

Cat. No.	Description	Qty
C2500	MacroTube®, 5mL, free of RNase & DNase, natural color, 4 bags of 50	200
C2500-AST	MacroTube®, 5mL, free of RNase & DNase, rainbow pack, 4 bags of 50	200
C2500-B	MacroTube®, 5mL, free of RNase & DNase, blue, 4 bags of 50	200
C2500-G	MacroTube®, 5mL, free of RNase & DNase, green, 4 bags of 50	200
C2500-OB	MacroTube®, 5mL, free of RNase & DNase, opaque black, 4 bags of 50	200
C2500-R	MacroTube®, 5mL, free of RNase & DNase, red, 4 bags of 50	200
C2500-Y	MacroTube®, 5mL, free of RNase & DNase, yellow, 4 bags of 50	200
C2520	MacroTube®, 5mL, sterile, free of RNase & DNase & pyrogens, 10 bags of 20	200
C2570	Adapter for 5mL tubes in 15mL centrifuge	4
C2581	Freezer box, polycarbonate, for 25 (5x5) 5mL tubes. (5.25 x 5.25 x 2.75in.)	5
C2586	CapLock Clip™, with break-away lifting tab	20
C2590	Two-tiered pipetting rack for 16 x 5mL MacroTubes®	1
C2592	High volume gripper rack for 50 x 5mL tubes	1
C2595	MyFuge™ 5 Mini-centrifuge for 4 x 5mL tubes, 5500rpm (see next page)	1



Adapter for 5mL tube  
in 15mL Centrifuge