

# BENCHWAVER™



- Undulating, 3-dimensional motion
- Adjustable speed and tilt angle
- Large, 13 inch platform
- Ideal for staining, washing & general mixing



B3D5000 with B3D5000-ST



B3D5000 with B3D5000-DIMP

The BenchWaver™ combines the basic rocking motion of traditional platform rockers with the circular motion of an orbital shaker. This combination is ideal for generating the optimum 3-dimensional "undulating" motion commonly preferred for molecular biology applications.

The tilt angle is instantly adjustable (no tools or disassembly required) from 0 to 10° and the mixing speed can be adjusted from 5 to 105 rpm, using the large, user-friendly LCD control panel. Combined control of both parameters (tilt and speed) makes the BenchWaver a versatile workhorse that can provide anything from gentle, slow moving waves to a rapid, more vigorous motion. Additionally, time of operation can be selected from 1 minute up to 99 hours or for continuous operation.

The large 13.5 x 13 inch platform is supplied complete with a non-slip, rubber mat. An optional stacking platform is available for users requiring additional mixing capacity. BenchWaver is a reliable and versatile instrument and can be used in refrigerated or incubated environments from 2° to 60°C.

## Technical Data:

Speed:	5 to 105 rpm
Tilt Angle:	Adjustable, 0° to 10°
Platform Dimension:	13.5 x 13 in. (34 x 33cm)
Dimensions:	13 x 13.5 x 9.75 in. 33 x 34.3 x 24.75 cm
Operating Temp. Range:	+2°C to +60°C
Load Capacity:	Up to 20lbs (9kg)
Weight:	5.5 kg / 12 lb
Electrical:	100-240V, 50-60 Hz, 50W
Warranty:	2 years

## Ordering Information:

B3D5000*	BenchWaver™ 3-D Rocker with flat mat platform, US plug
B3D5000-ST	Optional stacking platform (4.5" clearance) with stacking hardware
B3D5000-DIMP	Optional dimpled mat, for use with a variety of tubes

\* Includes US Plug. For EU plug, please add (-E)

**Benchmark**  
Scientific

**Benchmark Scientific Inc.**  
PO Box 709 Edison, NJ 08818  
Phone: 908-769-5555  
Fax: 908-222-1864  
Web: [www.BenchmarkScientific.com](http://www.BenchmarkScientific.com)  
Email: [Info@BenchmarkScientific.com](mailto:Info@BenchmarkScientific.com)