


# SureStand™

## MULTI-CHANNEL CAPABLE PIPETTE RACK



- **Universal - Fits all popular pipettes**
- **5mm acrylic - 40% thicker**
- **Non-skid rubber feet**
- **Multi-channel “window” slots**



SureStand's “window” slots are the first of their kind, allowing the user to store single channel or multi-channel pipettes easily and conveniently.

**MTC Bio**

2600A Main Street Extension  
Sayreville, NJ 08872  
Ph: 908-279-0280 www.mtcbio.com

P4403

P4406

3-Place Rack  
up to one multi-channel

SureStand™ Pipette racks are unique in their ability to accommodate multi-channel pipettes without sacrificing any of the stability of a traditional acrylic rack. The patent pending “windows” allow the user to insert the multi-channel head longways through the rack, while keeping the body of the pipette stable by securing the pipette at two points; under the finger hook, and near the bottom of the handle.

SureStand is available in 3 versions; a 3-place rack that accepts one multi-channel, a 5-place rack that accepts two multi-channels, and a 6-place rack that accepts four multichannels. Each SureStand is also outfitted with four non-skid rubber “feet” for extra stability. SureStand is guaranteed to be the sturdiest and most versatile pipette rack in your lab

P4405

5-Place Rack  
up to two multi-channels

### Ordering Information

Cat. No.	Description
P4403	SureStand Pipette Rack for 3 pipettes (up to one multi-channel)
P4405	SureStand Pipette Rack for 5 pipettes (up to two multi-channels)
P4406	SureStand Pipette Rack for 6 pipettes (up to four multi-channels)