

TECHNICAL DATA SHEET

Part Number	064035
Description	30 ml centrifuge tubes



Benefits :

- ✓ Tube material : polypropylene CoPolymer (PPC)
- ✓ Color : natural
- ✓ Cap material : HDPE
- ✓ Cap color : green
- ✓ Shape : freestanding
- ✓ Molded internally and feature molded-in graduations at 5 ml intervals
- ✓ Packed 20 tubes/tray

Certificate :

- ✓ Sterile, Gamma irradiated

Technical specifications :

- ✓ Capacity : 30 ml
- ✓ Withstand up to 6,000 x g
- ✓ Dimensions (Ø x h) : 30 x 88 mm

Sales unit :

480 Centrifuge tubes