

95-Pin Inoculator



95-Pin Inoculator

Our 95-Pin Inoculator is designed to pick up and transfer 95 samples into, or from, a regular 96-well format microplate or an inoculum tray.

The missing pin marks the H12 well location and serves as a blank control.

Inoculum Trays have the same footprint as regular microplate lids, and come packaged with the inoculator to be used in conjunction with it as a cover, or as a culture media plate. It can also be used as a general purpose microplate tray.

The 95-pin is ideal for antimicrobial susceptibility testing, urine cultures, and subcultures of bacterial colonies. It can also be used for cell/tissue cultures, cytotoxicity testing, serology dilution studies, and evaluation of growth media.

Material: Polystyrene

Gamma Sterilized

For more information, to request a sample or to contact your sales representative, please contact us:

800.421.6261

info@evergreensci.com

www.evergreensci.com

Catalog Number	Description	Packaging
222-8600-01A	95-Pin Inoculator Assembly	5/bag, 50/case
222-8023-01T	20mL Inoculum Tray	10/bag, 60/case
222-8021-01T	40mL Inoculum Tray	15/bag, 105/case

