

TECHNICAL DATA SHEET

Part Number	064227
Description	Specimen container



Benefits :

- ✓ Caplugs Evergreen wide-mouth container with screw caps
- ✓ Precision-molded from polypropylene
- ✓ Robust
- ✓ Lid color : white

Certificate :

- ✓ Non sterile
- ✓ Autoclavable
- ✓ The screwcaps are lined with XPE foam to help minimize evaporation during long-term storage and help eliminate leakage during shipping
- ✓ Chemically resistant to 95% ethyl alcohol and 10% formalin and some other commonly used organic solvents such as kerosene and isopropyl alcohol
- ✓ Can withstand boiling and freezing

Technical specifications :

- ✓ Volume : 950 ml
- ✓ Dimensions (Ø x h) : 120 x 138 mm

Sales unit :

24 Specimen containers