

### 24-Well Deep-Well Plates DP24 Series



- Sterile version is Nuclease (DNase and RNase) free
- Deep-Well plates have the same footprint as the regular microplates
- Well volume: 10 mL
- Flat (F) bottom with chamfered edges

### 48-Well Deep-Well Plates DP48 Series



- Sterile version is Nuclease (DNase and RNase) free
- Deep-Well plates have the same footprint as the regular microplates
- Well volume: 5 mL
- Flat (F) bottom with chamfered edges

### 96-Well Deep-Well Plates DP9, DP9C Series



- Sterile version is Nuclease (DNase and RNase) free
- Deep-Well plates have the same footprint as the regular microplates

## DP24 Series 24-Well Deep-Well Plates

Caplugs No.	Type	Material	Packaging	Carton Quantity
290-8573-08C	Standard	Polypropylene Co-Polymer (PPC)	5/bag	20
222-8685-08C	Sterile	Polypropylene Co-Polymer (PPC)	5/bag	20

### Lids for Plates

290-8022-03LS	Lids for Deep-Well Plates, 24-, 48- and 96-Well	Polypropylene (PP)	15/box	105
---------------	---	--------------------	--------	-----

## DP48 Series 48-Well Deep-Well Plates

Caplugs No.	Type	Material	Packaging	Carton Quantity
290-8560-08C	Standard	Polypropylene Co-Polymer (PPC)	5/bag	20
222-8672-08C	Sterile	Polypropylene Co-Polymer (PPC)	5/bag	20

### Lids for Plates

290-8022-03LS	Lids for Deep-Well Plates, 24-, 48- and 96-Well	Polypropylene (PP)	15/box	105
---------------	---	--------------------	--------	-----

## DP9, DP9C Series 96-Well Deep-Well Plates

Caplugs No.	Type	Well Volume	Bottom	Material	Packaging	Carton Quantity
290-8525-010	Standard	0.5 mL	Tapered (T)	Polystyrene (PS)	5/bag	20
222-8636-010	Sterile	0.5 mL	Tapered (T)	Polystyrene (PS)	5/bag	20
290-8538-080	Standard	0.5 mL	Tapered (T)	Polypropylene Co-Polymer (PPC)	5/bag	20
222-8649-080	Sterile	0.5 mL	Tapered (T)	Polypropylene Co-Polymer (PPC)	5/bag	20
290-8551-01N	Standard	1.0 mL	Tapered (T)	Polystyrene (PS)	5/bag	20
222-8662-01N	Sterile	1.0 mL	Tapered (T)	Polystyrene (PS)	5/bag	20
290-8553-08N	Standard	1.0 mL	Tapered (T)	Polypropylene Co-Polymer (PPC)	5/bag	20
222-8664-08N	Sterile	1.0 mL	Tapered (T)	Polypropylene Co-Polymer (PPC)	5/bag	20
290-8556-03C	Standard	2.0 mL	Flat (F) with chamfered edges	Polypropylene (PP)	5/bag	20
222-8668-08C	Sterile	2.0 mL	Flat (F) with chamfered edges	Polypropylene Co-Polymer (PPC)	5/bag	20

### Lids for Plates

290-8019-01L	Lids for 0.5 mL and 1.0 mL Deep-Well Plates	Polystyrene (PS)	15/box	105
290-8020-03L	Lids for 0.5 mL and 1.0 mL Deep-Well Plates	Polypropylene (PP)	15/bag	105