

# Slide-Sette Savers™



## Slide-Sette Savers™

For shipping and storing same patient pathology microscope slides and tissue/biopsy paraffin blocks in a single container. Quickly retrieve and send out samples for second opinions, additional examination, or long-term storage.

Slide-Sette Savers™ have a protective lid with living hinges. The inside is compartmented to accommodate both the slide and cassette in their own protective space.

Slide-Sette Saver™ is available in eight different colors to allow for easy organization by tissue type. Storage boxes, accommodating twenty savers each, are also available.

Catalog Number	Description	Packaging
240-3800-030	Slide-Sette Saver, natural	100/carton
240-3800-B30	Slide-Sette Saver, blue	100/carton
240-3800-K30	Slide-Sette Saver, fuchsia	100/carton
240-3800-R30	Slide-Sette Saver, red	100/carton
240-3800-P30	Slide-Sette Saver, pink	100/carton
240-3800-U30	Slide-Sette Saver, purple	100/carton
240-3800-Y30	Slide-Sette Saver, yellow	100/carton
240-3800-Z30	Slide-Sette Saver, black	100/carton

Catalog Number	Description	Packaging
240-3828-W80	Slide-Sette Saver Storage Box, white	3/pack
240-3829-Y80	Slide-Sette Saver Storage Box, yellow	3/pack



For more information, or to order samples, please call or email:



[www.evergreensci.com](http://www.evergreensci.com)

800-421-6261

[info@evergreensci.com](mailto:info@evergreensci.com)