

Slide-Sette Savers™



Slide-Sette Savers™

For shipping and storing same patient pathology microscope slides and tissue/biopsy paraffin blocks in a single container. Quickly retrieve and send out samples for second opinions, additional examination, or long-term storage.

Slide-Sette Savers™ have a protective lid with living hinges. The inside is compartmented to accommodate both the slide and cassette in their own protective space.

Slide-Sette Saver™ is available in eight different colors to allow for easy organization by tissue type. Storage boxes, accommodating twenty savers each, are also available.

| Catalog Number | Description | Packaging |
|----------------|----------------------------|------------|
| 240-3800-030 | Slide-Sette Saver, natural | 100/carton |
| 240-3800-B30 | Slide-Sette Saver, blue | 100/carton |
| 240-3800-K30 | Slide-Sette Saver, fuchsia | 100/carton |
| 240-3800-R30 | Slide-Sette Saver, red | 100/carton |
| 240-3800-P30 | Slide-Sette Saver, pink | 100/carton |
| 240-3800-U30 | Slide-Sette Saver, purple | 100/carton |
| 240-3800-Y30 | Slide-Sette Saver, yellow | 100/carton |
| 240-3800-Z30 | Slide-Sette Saver, black | 100/carton |

| Catalog Number | Description | Packaging |
|----------------|---------------------------------------|-----------|
| 240-3828-W80 | Slide-Sette Saver Storage Box, white | 3/pack |
| 240-3829-Y80 | Slide-Sette Saver Storage Box, yellow | 3/pack |



For more information, or to order samples, please call or email:



www.evergreensci.com

800-421-6261

info@evergreensci.com