

KIMBLE® KONTES® Sterile Closed-System Tissue Grinder, with Molded Tip, 50 mL



KIMBLE® KONTES® Sterile Closed-System Tissue Grinder, with Molded Tip, 50 mL

- Manufactured from polypropylene (tubes); pestles are 30% glass-filled PP and caps are HDPE
- Sterile Tissue Grinder available with two different tips: molded (closed system 1) and vitrified (closed system 2)
- The proprietary silica casting process in Closed System 2 tissue grinders results in a glass-like abrasive tip surface that makes it easier to homogenize fibrous samples
- Tissues are ground within a sealed container to minimize the risk of personal contact -

exposure to sample aerosols is eliminated

- The molded-in abrasive surface on the pestle tip leaves no sediment to obstruct sample examination
- Narrow pestles eliminate the potential of sample overflow common with other types of grinders
- All are gamma-sterilized in easy-open Tyvek packs
- Each grinder comes with a fully assembled pestle, a conical sample tube with solid top cap, an adhesive identification label and an instruction card.

Catalog Number: **749600-0050**

DOMINIQUE DUTSCHER SAS