


# Pipette Accessories



*Gilson pipette accessories  
are the perfect addition  
to your pipette collection*

- > **To ensure safety and convenience as well as ISO compliance**, pipettes must be protected from possible damage, contamination or chemical attack and stored safely in an upright position.
- > In laboratories, where **GLP and sample integrity are equally important**, you need solutions that work for your pipetting system.
- > Besides being convenient, Gilson accessories provide **solutions for storage, handling, tracing and identification** to make your path to discovery easier.



# Accessories

## ■ Carrousel™ Pipette Stand

This rotary carrousel rack holds up to seven Gilson Pipettes, all models with tips fitted except Pipetman Concept®. Carrousel is robust, stable and space saving, making your work bench look neat and tidy. Its heavy base prevents toppling over.



## ■ Trio™ Pipette Stand

The Trio stand features enhanced quality and reliability. Its angled, transparent design provides a convenient and secure way to store your Gilson pipettes. The transparent body allows you to see the objects placed behind.



## ■ Single™ Pipette Holder

Single is designed to conveniently hold all Gilson Pipettes except Pipetman Concept where space is limited (fume hoods, sterile hoods, etc.). Stuck to any smooth, dry surface or fastened with screws, it keeps your pipette within reach. (Chemicals and UV light resistant).



## ■ Pipetman Concept® Rapid-Charge Stand

This robust stand stores the Pipetman Concept in an upright position while charging 80% in one hour, fully in two. A second battery block can also be kept charged and available for continuous operation.



## ■ Gilson Reagent Reservoirs

Disposable reagent reservoirs are ideal for dispensing reagents with multichannel pipettes. The Gilson Reagent reservoirs are available in two sizes, 25 mL and 50 mL, standard or sterilized. They meet stringent criteria of reliability and quality.



- > Wide, stable base
- > Angled sides and V-shaped bottoms
- > Pour-off spouts on all corners
- > Graduations on inside walls at the middle volume
- > Made of white polystyrene
- > Sterilized via Gamma irradiation and individually wrapped

## ■ Coloris™ Colored Identification Clips for Pipetman®

GLP-friendly, these brightly colored clips are used to identify or personalize your Pipetman pipettes. Just write your name or application to avoid mix-ups. Made of polyacetal, Coloris Clips accept permanent ink and just snap into place.



## ■ Pipetman® Ultra & Concept® Name Tags

Personalize your Pipetman Ultra and Pipetman Concept by writing your name or application on one of the colored stickers (protected by the plastic window).



### Description

Description	Part Number
Carrousel for 7 pipettes	F161401
Trio, stand for 3 pipettes	F161405
Single pipette holder	F161406
Rapid-Charge Stand, PIPETMAN CONCEPT	F30750
Reagent Reservoir 25ml, 10 packs of 10 reservoirs	F267660
Reagent Reservoir 50ml, 10 packs of 10 reservoirs	F267670
Sterilized Reagent Reservoir 25ml, 50 individually wrapped	F267661
Sterilized Reagent Reservoir 50ml, 50 individually wrapped	F267671
10 Coloris clips, mixed colours, PIPETMAN	F161301
10 Coloris clips, red, PIPETMAN	F161302
10 Coloris clips, yellow, PIPETMAN	F161303
10 Coloris clips, green, PIPETMAN	F161304
10 Coloris clips, blue, PIPETMAN	F161305
10 Coloris clips, white, PIPETMAN	F161306
Set of 5 name tags, PIPETMAN ULTRA	F2070901
Set of 5 name tags, PIPETMAN CONCEPT	F3071107

Distributed by

[www.gilson.com](http://www.gilson.com) | [sales@gilson.com](mailto:sales@gilson.com) | [service@gilson.com](mailto:service@gilson.com) | [training@gilson.com](mailto:training@gilson.com)

**Gilson, Inc. World Headquarters** | 3000 Parmentier Street | P.O. Box 620027 | Middleton, WI 53562-0027 USA  
Tel: (1) 800-445-7661 or (1) 608-836-1551 | Fax: (1) 608-831-4451

**Gilson S.A.S.** | 19, avenue des Entrepreneurs | BP 145, F-95400 Villiers-le-Bel, FRANCE  
Tel: (33-1) 34 29 50 00 | Fax: (33-1) 34 29 50 20

**ISO9001 Certified**

L1800479E - © 2006 Gilson S.A.S. - All rights reserved - Printed in France, June 2006.  
Specifications subject to change without notice. Gilson, Inc. and Gilson S.A.S are ISO 9001 certified.

