

distriman®

FASTER MULTIPLE DISPENSING

DISTRIMAN® is a positive-displacement, continuously adjustable, repetitive pipette with direct read out designed to simplify multiple dispensing when used with DISTRITIPS syringes. With DISTRIMAN, you can dispense the exact volume required for any protocol, even fractional volumes. It is ideal for repetitive dispensing in clinical, hospital, biological, chemical, food and beverage, pharmaceutical, cosmetics, and forensics labs.

SIMPLIFIES MULTIPLE DISPENSING

Dispense exact volumes without tedious calculations and ensure optimal pipetting results even with volatile or viscous liquids.

SAVE MONEY

DISTRIMAN covers a large volume range with only three DISTRITIPS syringes, making it very economical. DISTRITIPS are available in standard or sterilized form (boxes of 50).

DISTRIMAN is the only continuously adjustable mechanical repetitive pipette with direct volume readout available on the market.



DISTRIMAN®

Description	Part Number
DISTRIMAN (1 µL to 1.25 mL)	F164001

DISTRITIPS

DISTRITIPS (Pack of 50)	Aliquot Range	Part Number	Volume (µL)	Maximum Permissible Errors			
				Gilson		ISO 8655	
				Systematic Error (µL)	Random Error (µL)	Systematic Error (µL)	Random Error (µL)
125 µL Micro Micro ST*	1-12.5 µL	F164100 F164130	2	± 0.100	≤ 0.080	± 0.2	≤ 0.1
			5	± 0.125	≤ 0.075	± 0.2	≤ 0.1
			10	± 0.200	≤ 0.100	± 0.2	≤ 0.1
1250 µL Mini Mini ST*	10-125 µL	F164110 F164140	20	± 0.80	≤ 0.20	± 1.5	≤ 0.6
			50	± 1.00	≤ 0.40	± 1.5	≤ 0.6
			100	± 1.00	≤ 0.60	± 1.5	≤ 0.6
12.5 mL Maxi Maxi ST*	0.1-1.25 mL	F164120 F164150	200	± 6.0	≤ 1.0	± 12.0	≤ 4.0
			500	± 7.5	≤ 1.5	± 12.0	≤ 4.0
			1000	± 10.0	≤ 2.5	± 12.0	≤ 4.0

*ST means sterilized syringes (individually wrapped)