

Declaration of Conformity

Manufacturer: **Gilson SAS** - 19 avenue des Entrepreneurs - BP 145 - 95400 Villiers-le-Bel - France

Type of Equipment: **Laboratory Equipment**


Application of Council Directives: **89/336/EEC, 73/23/EEC and 2002/96/EEC**

Standards to which Conformity is Declared:

Products		EN61010-1 UL61010-1 CAN/CSA-C22.2 N°61010-1	EN50081-1 EN50082-1	EN50081-2 EN50082-2	EN61326
Autoinjector:	234	x	x	x	x
High Pressure Pumps:	305, 306, 307, 321, 322, 331, 332, 333, 334	x x	x x	x x	x x
Peristaltic Pump:	Minipuls 3 Minipuls Evolution	x x	x	x	x x
Dilutors:	402 Syringe Pump, 402 Dilutor Dispenser,	x x	x x	x x	x x
Manometric Modules:	805, 806, 807	x	x	x	x
Temperature Regulators:	832	x	x	x	x

I, the undersigned, hereby declare that the equipment specified above, conforms to the above Directives and Standards.

Hervé Le Dorze
Quality Assurance Manager
Place: Villiers le Bel
Issue date: January 25, 2007



LT801354/K

The Gilson SAS quality assurance system is registered as complying with the ISO 9001 International Quality Standard. Registration Certificates N°6621 and N°18623 BVQI.