



CSP Dairy™ Single Pack Dipper

Single Use Clean Dipper for Dairy Sampling

Applications

Fields of use

- Field collection of milk and related dairy products
- Payment and quality control

Options

- Standard dipper is colourless – colours upon request

Benefits

Clean Single Use Dipper

- Individually packed in sealed perforated bag
- Translucent in colour for visual inspection
- Pure Polypropylene suitable for food contact
- 100% Recyclable after use

Regulatory and Quality

- ISO-9001 certified company
- Complies with Regulation (EC) No 1935/2004

Characteristics

Please note that these values are average measurements
Tolerance +/- 0.5 mm (will vary based on product)

Internal volume	30 ml
Overall length	200 mm
Dipper handle length	164,2 mm
Dipper handle width	13,8 mm
Body height	40,8 mm
Top body diameter	52,1 mm
Bottom body diameter	28,3 mm
Weight	8,6 grams

Production mould # 706

Packaging

ART. # 03 05 2200 MD00 CN

Clear in color
370 Pieces per box
30 Boxes
11.100 Pieces per pallet

